

# LIBRARY PREPARATION

## NEBNext<sup>®</sup> Multiplex Oligos for Illumina<sup>®</sup> (96 Index Primers)

Instruction Manual

NEB #E6609S/L  
96/384 reactions  
Version 1.3 6/16






USER™ is protected by U.S. Patent No. 7,435,572 (New England Biolabs, Inc.).

NEW ENGLAND BIOLABS® and NEBNEXT® are registered trademarks of New England Biolabs, Inc. USER™ is a trademark of New England Biolabs, Inc.

AMPURE® is a registered trademark of Beckman Coulter, Inc.  
BIOANALYZER® is a registered trademark of Agilent Technologies, Inc.  
EPPENDORF® and LOBIND® are registered trademarks of Eppendorf AG.  
ILLUMINA®, HISEQ® and MISEQ® are registered trademarks of Illumina, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

© Copyright 2016, New England Biolabs, Inc; all rights reserved.



## Table of Contents:

---

Applications.....	2
Chapter 1	
Setting up the PCR Reaction.....	5
Chapter 2	
Index Pooling Guidelines.....	8
NEBNext Adaptor for Illumina.....	14
USER™ Enzyme.....	15
NEBNext Index/Universal Primer Mix Plate.....	16-21
Revision History.....	22

## The NEBNext Multiplex Oligos for Illumina (96 Index Primers) Includes:

---

*The volumes provided are sufficient for preparation of up to 96 reactions (NEB #E6609S) and 384 reactions (NEB #E6609L). All reagents should be stored at -20°C.*

NEBNext Adaptor for Illumina

USER Enzyme

NEBNext Index/Universal Primer Mix Plate

Each well contains the Universal PCR Primer plus one of the Index Primers

## Required Materials Not Included:

---

Enzymes and buffers appropriate for DNA, ChIP or mRNA Library Preparation (NEB #E7645, E7595, E7420, E7530, E7370, E7445, E6040, E6110, E6240, E6056, E6000, E6200, E6100)

SPRIselect or AMPure® XP Beads (Beckman Coulter, Inc.)

DNA LoBind® Tubes (Eppendorf®)/PCR Plates

Magnetic Stand

Bioanalyzer® (Agilent Technologies, Inc.)

Freshly Prepared 80% Ethanol

## Applications:

The NEBNext Multiplex Oligos for Illumina (96 Index Primers) contains adaptors and primers that are ideally suited for multiplex sample preparation for next-generation sequencing on the Illumina platform (Illumina, Inc.). Each of these components must pass rigorous quality control standards and is lot controlled, both individually and as a set of reagents.

**Lot Control:** The lots provided in the NEBNext Multiplex Oligos for Illumina (96 Index Primers) are managed separately and are qualified by additional functional validation. Individual reagents undergo standard enzyme activity and quality control assays, and also meet stringent criteria in the additional quality controls listed on each individual component page.


**Functionally Validated:** Each set of reagents is functionally validated together through construction and sequencing of genomic DNA libraries on the Illumina platform.

For larger volume requirements, customized and bulk packaging is available by purchasing through the OEM/Bulks department at NEB. Please contact [OEM@neb.com](mailto:OEM@neb.com) for further information.

## Workflow Overview:

Designed for use in library prep for DNA, ChIP DNA and RNA (but not Small RNA), the NEBNext Adaptors enable high-efficiency adaptor ligation and high library yields, with minimized adaptor-dimer formation. Incorporating a novel hairpin loop structure, the NEBNext Adaptor ligates with increased efficiency to end-repaired, dA-tailed DNA. The loop contains a U, which is removed by treatment with USER Enzyme (a combination of UDG and Endo VIII), to open up the loop and make it available as a substrate for PCR. During PCR, barcodes can be incorporated by use of the NEBNext index primers, thereby enabling multiplexing. The 96 8-base index primers included in this kit are pre-mixed with the universal primer and are packaged in a single-use 96-well plate with a pierceable foil seal. NEBNext Oligos can be used with NEBNext products, and with other standard Illumina-compatible library preparation protocols.

Figure 1: Workflow demonstrating the use of NEBNext Multiplex Oligos for Illumina (96 Index Primers)



## Please Refer to the Kit Specific Protocol for using the NEBNext Multiplex Oligos for Illumina:

The following kits are designed for use with the NEBNext Multiplex Oligos for Illumina:

- #E7645, NEBNext Ultra II DNA Library Prep Kit for Illumina
- #E7595, NEBNext Ultra II Ligation Module
- #E7420, NEBNext Ultra Directional RNA Library Prep Kit for Illumina
- #E7530, NEBNext Ultra RNA Library Prep Kit for Illumina
- #E7370, NEBNext Ultra DNA Library Prep Kit for Illumina
- #E7445, NEBNext Ultra Ligation Module
- #E6040, NEBNext DNA Library Prep Master Mix Set for Illumina
- #E6110, NEBNext mRNA Library Prep Master Mix Set for Illumina
- #E6240, NEBNext ChIP-Seq Library Prep Master Mix Set for Illumina
- #E6056, NEBNext Quick Ligation Module
- #E6000, NEBNext DNA Library Prep Reagent Set for Illumina
- #E6200, NEBNext ChIP-Seq Library Prep Reagent Set for Illumina
- #E6100, NEBNext mRNA Library Prep Reagent Set for Illumina


# 1

## Setting up the PCR Reaction

---

NEBNext Multiplex Oligos for Illumina  
(96 Index Primers)  
Instruction Manual

## Symbols

 This caution sign signifies a step in the protocol that has multiple paths leading to the same end point but is dependent on a user variable, like the amount of input DNA.

### 1.1. PCR Amplification


 For < 96 samples, follow the protocol in Section 1.1A. For 96 samples, follow the protocol in Section 1.1B.

#### 1.1A. Setting up the PCR reactions (< 96 samples)

- 1.1A.1. Determine the number of libraries that will be amplified and pooled for subsequent sequencing.
- 1.1A.2. Ensure that you choose a valid combination of barcode primers based on color balance guidelines in Chapter 2.
- 1.1A.3. Thaw the NEBNext Index/Universal Primer Mix plate for 10-15 minutes at room temperature.
- 1.1A.4. Remove the hard plastic plate cover. Briefly centrifuge the plate (280 ×g for ~1 min) to collect all of the primer at the bottom of each well.
- 1.1A.5. Orient the NEBNext Index/Universal Primer Mix Plate as indicated in Figure 1.1 (red stripe towards the user). With a pipette tip, pierce the desired well(s) (Figure 1.1A) and transfer the volume of primer mix required for the PCR reaction to the PCR plate/tubes (see specific library construction manual for protocol). It is important to change pipette tips before piercing a new well to avoid cross contamination of indexed primers. Alternatively, the wells can be pierced using the bottom of clean PCR strip tubes (see Figure 1.1B) prior to pipetting the primer mix. Use a new, clean strip tube for each new well to be pierced.  
*Note: Each well contains the Universal Primer and the Index Primer. There is enough primer in each well for one PCR reaction. Do not reuse primer if the seal has been previously pierced to avoid contamination with other indexed primers.*
- 1.1A.6. Proceed with the PCR reaction according to the specific library construction manual.



Figure 1.1: NEBNext Index/Universal Primer Mix plate.



### 1.1B. Setting up the PCR reactions (96 samples)

- 1.1B.1. Thaw the NEBNext Index/Universal Primer Mix plate for 10-15 minutes at room temperature.
- 1.1B.2. Remove the hard plastic plate cover. Briefly centrifuge the plate (280  $\times$ g for ~1 min) to collect all of the primer at the bottom of each well.
- 1.1B.3. Orient the NEBNext Index/Universal Primer Plate as indicated in Figure 1.1 (red stripe towards the user). With a pipette tip, pierce the wells (Figure 1.1A) and transfer the volume of primer mix required for the PCR reaction to the PCR plate (see specific library construction manual for protocol). It is important to change pipette tips before piercing a new well to avoid cross contamination of indexed primers. Alternatively, the wells can be pierced using the bottom of clean PCR strip tubes (see Figure 1.1B) prior to pipetting the primer mix. Use a new, clean strip tube for each new well to be pierced.

*Note: Each well contains the Universal Primer and the Index Primer. There is enough primer in each well for one PCR reaction. Do not reuse primer if the seal has been previously pierced to avoid contamination with other indexed primers.*

- 1.1B.4. Proceed with the PCR reaction according to the specific library construction manual.

# 2

## Index Pooling Guidelines

---

NEBNext Multiplex Oligos for Illumina  
(96 Index Primers)  
Instruction Manual

**For a link to download a sample sheet with the index sequences for use with the Illumina Experiment Manager (IEM) please go to our FAQ's tab on [www.neb.com/E6609](http://www.neb.com/E6609) – NEBNext Multiplex Oligos for Illumina (96 Index Primers) (NEB #E6609).**

For the HiSeq®/MiSeq®, Illumina uses a red laser/LED to sequence bases A and C and a green laser/LED to sequence bases G and T. For each cycle, both the red and the green channel need to be read to ensure proper image registration (i.e. A or C must be in each cycle, and G or T must be in each cycle). If this color balance is not maintained, sequencing the index read could fail. Table 2.1 lists some valid combinations (up to 4-plex) that can be sequenced together. For combinations > 4 choose any 4-plex combination and add any other index as needed.

Table 2.1

PLEX	INDEX PRIMER
2	P18-B6 and P22-B10 P37-D1 and P42-D6 P52-E4 and P76-G4 P68-F8 and P95-H11
3	P4-A4, P11-A11, and P12-A12 P15-B3, P22-B10, and, P24-B12 P25-C1, P31-C7, and P33-C9 P37-D1, P42-D6, and, P48-D12 P49-E1, P54-E6, and, P55-E7 P64-F4, P69-F9, and P71-F11 P76-G4, P77-G5, and P83-G11 P87-H3, P93-H9, and P94-H10

Table 2.1 (continued)

PLEX	INDEX PRIMER
4	P1-A1, P2-A2, P3-A3, and P4-A4
	P5-A5, P6-A6, P8-A8, and P10-A10
	P13-B1, P14-B2, P15-B3, and P16-B4
	P17-B5, P18-B6, P19-B7, and P20-B8
	P25-C1, P26-C2, P27-C3, and P30-C6
	P28-C4, P29-C5, P32-C8, and P35-C11
	P37-D1, P38-D2, P39-D3, and P40-D4
	P45-D9, P46-D10, P47-D11, and P48-D12
	P49-E1, P50-E2, P51-E3, and P52-E4
	P56-E8, P58-E10, P59-E11, and P60-E12
	P61-F1, P62-F2, P63-F3, and P69-F9
	P64-F4, P65-F5, P66-F6, and P67-F7
	P73-G1, P74-G2, P75-G3, and P76-G4
	P80-G8, P82-G10, P83-G11, and P84-G12
	P85-H1, P86-H2, P87-H3, and P89-H5
	P91-H7, P94-H10, P95-H11, and P96-H12

Table 2.2 lists each index sequence color coded to correspond to the red/green channel. For combinations of valid indices, ensure that you will have signal in both the red and green channels in each cycle. See examples below:

GOOD																	
PRIMER	INDEX SEQUENCE							PRIMER	INDEX SEQUENCE								
P1-A1	T	T	A	C	C	G	A	C	P41-D5	G	A	C	G	T	C	A	T
P2-A2	A	G	T	G	A	C	C	T	P42-D6	C	T	T	A	C	A	G	C
P3-A3	T	C	G	G	A	T	T	C	P43-D7	T	C	C	A	T	T	G	C
P4-A4	C	A	A	G	G	T	A	C	P44-D8	A	G	C	G	A	G	A	T
	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓

BAD																	
PRIMER	INDEX SEQUENCE							PRIMER	INDEX SEQUENCE								
P9-A9	C	G	C	A	A	C	T	A	P56-E8	T	A	T	G	G	C	A	C
P10-A10	C	G	T	A	T	C	T	C	P57-E9	C	T	C	G	A	A	C	A
P11-A11	G	T	A	C	A	C	C	T	P58-E10	C	A	A	C	T	C	C	A
P12-A12	C	G	G	C	A	T	T	A	P59-E11	G	T	C	A	T	C	G	T
	✓	✗	✓	✗	✓	✓	✓	✓		✓	✓	✓	✓	✓	✗	✓	✓

Table 2.2 Index Sequences

INDEX PRIMER	INDEX SEQUENCE READ							
P1-A1	T	T	A	C	C	G	A	C
P2-A2	A	G	T	G	A	C	C	T
P3-A3	T	C	G	G	A	T	T	C
P4-A4	C	A	A	G	G	T	A	C
P5-A5	T	C	C	T	C	A	T	G
P6-A6	G	T	C	A	G	T	C	A
P7-A7	C	G	A	A	T	A	C	G
P8-A8	T	C	T	A	G	G	A	G
P9-A9	C	G	C	A	A	C	T	A
P10-A10	C	G	T	A	T	C	T	C
P11-A11	G	T	A	C	A	C	C	T
P12-A12	C	G	G	C	A	T	T	A
P13-B1	T	C	G	T	C	T	G	A
P14-B2	A	G	C	C	T	A	T	C
P15-B3	C	T	G	T	A	C	C	A
P16-B4	A	G	A	C	C	T	T	G
P17-B5	A	G	G	A	T	A	G	C
P18-B6	C	C	T	T	C	C	A	T
P19-B7	G	T	C	C	T	T	G	A
P20-B8	T	G	C	G	T	A	A	C
P21-B9	C	A	C	A	G	A	C	T
P22-B10	T	T	A	C	G	T	G	C
P23-B11	C	C	A	A	G	G	T	T
P24-B12	C	A	C	G	C	A	A	T
P25-C1	T	T	C	C	A	G	G	T
P26-C2	T	C	A	T	C	T	C	C
P27-C3	G	A	G	A	G	T	A	C
P28-C4	G	T	C	G	T	T	A	C
P29-C5	G	G	A	G	G	A	A	T
P30-C6	A	G	G	A	A	C	A	C
P31-C7	C	A	G	T	G	C	T	T
P32-C8	C	T	T	G	C	T	A	G

INDEX PRIMER	INDEX SEQUENCE READ							
P33-C9	T	G	G	A	A	G	C	A
P34-C10	A	G	C	T	A	A	G	C
P35-C11	G	A	A	C	G	G	T	T
P36-C12	G	G	A	A	T	G	T	C
P37-D1	T	A	C	G	G	T	C	T
P38-D2	C	C	A	G	T	A	T	C
P39-D3	T	C	T	A	C	G	C	A
P40-D4	G	T	A	A	C	C	G	A
P41-D5	G	A	C	G	T	C	A	T
P42-D6	C	T	T	A	C	A	G	C
P43-D7	T	C	C	A	T	T	G	C
P44-D8	A	G	C	G	A	G	A	T
P45-D9	C	A	A	T	A	G	C	C
P46-D10	A	A	G	A	C	A	C	C
P47-D11	C	C	A	G	T	T	G	A
P48-D12	T	G	G	T	G	A	A	G
P49-E1	A	A	G	A	C	C	G	T
P50-E2	T	T	G	C	G	A	G	A
P51-E3	G	C	A	A	T	T	C	C
P52-E4	G	A	A	T	C	C	G	T
P53-E5	C	C	G	C	T	T	A	A
P54-E6	T	A	C	C	T	G	C	A
P55-E7	G	T	C	G	A	T	T	G
P56-E8	T	A	T	G	G	C	A	C
P57-E9	C	T	C	G	A	A	C	A
P58-E10	C	A	A	C	T	C	C	A
P59-E11	G	T	C	A	T	C	G	T
P60-E12	G	G	A	C	A	T	C	A
P61-F1	C	A	G	G	T	T	C	A
P62-F2	G	A	A	C	G	A	A	G
P63-F3	C	T	C	A	G	A	A	G
P64-F4	C	A	T	G	A	G	C	A

INDEX PRIMER	INDEX SEQUENCE READ							
P65-F5	G	A	C	G	A	A	C	T
P66-F6	A	G	A	C	G	C	T	A
P67-F7	A	T	A	A	C	G	C	C
P68-F8	G	A	A	T	C	A	C	C
P69-F9	G	G	C	A	A	G	T	T
P70-F10	G	A	T	C	T	T	G	C
P71-F11	C	A	A	T	G	C	G	A
P72-F12	G	G	T	G	T	A	C	A
P73-G1	T	A	G	G	A	G	C	T
P74-G2	C	G	A	A	T	T	G	C
P75-G3	G	T	C	C	T	A	A	G
P76-G4	C	T	T	A	G	G	A	C
P77-G5	T	C	C	A	C	G	T	T
P78-G6	C	A	A	C	A	C	A	G
P79-G7	G	C	C	T	T	A	A	C
P80-G8	G	T	A	A	G	G	T	G
P81-G9	A	G	C	T	A	C	C	A
P82-G10	C	T	T	C	A	C	T	G
P83-G11	G	G	T	T	G	A	A	C
P84-G12	G	A	T	A	G	C	C	A
P85-H1	T	A	C	T	C	C	A	G
P86-H2	G	G	A	A	G	A	G	A
P87-H3	G	C	G	T	T	A	G	A
P88-H4	A	T	C	T	G	A	C	C
P89-H5	A	A	C	C	A	G	A	G
P90-H6	G	T	A	C	C	A	C	A
P91-H7	G	G	T	A	T	A	G	G
P92-H8	C	G	A	G	A	G	A	A
P93-H9	C	A	G	C	A	T	A	C
P94-H10	C	T	C	G	A	C	T	T
P95-H11	C	T	T	C	G	G	T	T
P96-H12	C	C	A	C	A	A	C	A

# NEBNext Adaptor for Illumina

**#E6612A: 0.96 ml**

**Concentration: 15  $\mu$ M**

Note: Large kit contains 4 x 0.96 ml

5'-/5Phos/GAT CGG AAG AGC ACA CGT CTG AAC TCC AGT C/ideoxyU/A CAC TCT TTC CCT ACA CGA  
CGC TCT TCC GAT C-s-T-3'

Where -s- indicates phosphorothioate bond.

**Store at -20°C**

## Quality Control Assays

**16-Hour Incubation:** 50  $\mu$ l reactions containing this adaptor and 1  $\mu$ g of HindIII digested Lambda DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50  $\mu$ l reactions containing this reaction buffer at 1X concentration and 1  $\mu$ g T3 DNA incubated for 16 hours at 37°C also results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

**Endonuclease Activity:** Incubation of a minimum of 5  $\mu$ l of this adaptor with 1  $\mu$ g of  $\phi$ X174 RF 1 DNA in assay buffer for 4 hours at 37°C in 50  $\mu$ l reactions results in < 10% conversion to RF II as determined by agarose gel electrophoresis.

**Phosphatase Activity:** Incubation of a minimum of 10  $\mu$ l of this adaptor in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl<sub>2</sub>) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**RNase Activity:** Incubation of this adaptor with 40 ng of a FAM-labeled RNA transcript for 16 hours at 37°C results in no detectable RNase activity as determined by polyacrylamide gel electrophoresis.

**Lot Controlled**



## USER Enzyme

#E6610A: 0.288 ml

#E6610AA: 0.576 ml

(2 vials)

Store at  $-20^{\circ}\text{C}$

Supplied in: 50 mM KCl, 5 mM NaCl, 10 mM Tris-HCl (pH 7.4 @  $25^{\circ}\text{C}$ ), 0.1 mM EDTA, 1 mM DTT, 175  $\mu\text{g/ml}$  BSA and 50% Glycerol

### Quality Control Assays

**Non-Specific DNase Activity (16 Hour):** A 50  $\mu\text{l}$  reaction in NEBuffer 1 containing 1  $\mu\text{g}$  of Lambda DNA and a minimum of 50 units of Uracil DNA Glycosylase incubated for 16 hours at  $37^{\circ}\text{C}$  results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis. A 50  $\mu\text{l}$  reaction in Endonuclease VIII Reaction Buffer containing 1  $\mu\text{g}$  of Lambda-HindIII DNA and a minimum of 25 units of Endonuclease VIII incubated for 16 hours at  $37^{\circ}\text{C}$  results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

**Exonuclease Activity (Radioactivity Release):** A 50  $\mu\text{l}$  reaction in NEBuffer 1 containing 1  $\mu\text{g}$  of a mixture of single and double-stranded [ $^3\text{H}$ ] *E. coli* DNA and a minimum of 50 units of Uracil DNA Glycosylase incubated for 4 hours at  $37^{\circ}\text{C}$  releases < 0.1% of the total radioactivity. A 50  $\mu\text{l}$  reaction in Endonuclease VIII Reaction Buffer containing 1  $\mu\text{g}$  of a mixture of single and double-stranded [ $^3\text{H}$ ] *E. coli* DNA and a minimum of 10 units of Endonuclease VIII incubated for 4 hours at  $37^{\circ}\text{C}$  releases < 0.5% of the total radioactivity.

**Endonuclease Activity (Nicking):** A 50  $\mu\text{l}$  reaction in UDG Reaction Buffer containing 1  $\mu\text{g}$  of supercoiled  $\phi\text{X174}$  DNA and a minimum of 50 units of Uracil DNA Glycosylase incubated for 4 hours at  $37^{\circ}\text{C}$  results in < 10% conversion to the nicked form as determined by agarose gel electrophoresis.

**Phosphatase Activity:** Incubation of a minimum of 10  $\mu\text{l}$  of USER at a 1X concentration in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM  $\text{MgCl}_2$ ) containing 2.5 mM *p*-nitrophenyl phosphate at  $37^{\circ}\text{C}$  for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**Lot Controlled**

# NEBNext Index/Universal Primer Mix

#E6611A: 1 plate (10 µl/well)

Concentration: 5 µM each primer

Note: Large kit contains 4 plates (10 µl/well).

**Description:** 96 Index Primers & the Universal Primer are included for producing barcoded libraries.

NEBNext Universal PCR Primer for Illumina:

5'-AAT GAT ACG GCG ACC ACC GAG ATC TAC ACT CTT TCC CTA CAC GAC GCT CTT CCG ATC-s-T-3'

Store at -20°C

Where -s- indicates phosphorothioate bond.

INDEX PRIMER	INDEX PRIMER SEQUENCE	EXPECTED INDEX SEQUENCE READ
P1-A1	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GTCGGTAA</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TTACCGAC
P2-A2	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AGGTCACT</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	AGTGACCT
P3-A3	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GAATCCGAG</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TCGGATTC
P4-A4	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GTACCTTG</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CAAGGTAC
P5-A5	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CATGAGGAG</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TCCTCATG
P6-A6	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGACTGAC</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GTCAGTCA
P7-A7	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CGTATTCG</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CGAATACG
P8-A8	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CTCCTAGAG</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TCTAGGAG
P9-A9	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TAGTTGCG</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CGCAACTA
P10-A10	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GAGATACG</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CGTATCTC
P11-A11	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AGGTGTAC</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GTACACCT
P12-A12	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TAATGCCG</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CGGCATTA
P13-B1	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TCAGACGA</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TCGTCTGA
P14-B2	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GATAGGCT</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	AGCCTATC
P15-B3	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGGTACAG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CTGTACCA
P16-B4	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CAAGGCTC</b> TGTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	AGACCTTG
P17-B5	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GCTATCCT</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	AGGATAGC
P18-B6	5'-CAAGCAGAAGACGGCATAACGAGAT <b>ATGGAGGG</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CCTTCCAT

## NEBNext Index/Universal Primer Mix (Cont.)

INDEX PRIMER	INDEX PRIMER SEQUENCE	EXPECTED INDEX SEQUENCE READ
P19-B7	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TC AAGGAC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GTCCTTGA
P20-B8	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GTTACGCA</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TGCGTAAC
P21-B9	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AGTCTGTG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CACAGACT
P22-B10	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GCACGTAAG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TTACGTGC
P23-B11	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AACCTTGGG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CCAAGGTT
P24-B12	5'-CAAGCAGAAGACGGCATAACGAGAT <b>ATTGCGTG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CACGCAAT
P25-C1	5'-CAAGCAGAAGACGGCATAACGAGAT <b>ACCTGGAAG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TTCCAGGT
P26-C2	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GGAGATGAG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TCATCTCC
P27-C3	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GTACTCTC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GAGAGTAC
P28-C4	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GTAACGAC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GTCGTTAC
P29-C5	5'-CAAGCAGAAGACGGCATAACGAGAT <b>ATTCCTCCG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GGAGGAAT
P30-C6	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GTGTTCTC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	AGGAACAC
P31-C7	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AAGCACTG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CAGTGCTT
P32-C8	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CTAGCAAG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CTTGCTAG
P33-C9	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGCTTCCA</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TGGAAGCA
P34-C10	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GCTTAGCT</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	AGCTAAGC
P35-C11	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AACCGTTC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GAACGGTT
P36-C12	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GACATTCCG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GGAATGTC
P37-D1	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AGACCGTA</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TACGGTCT
P38-D2	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GATACTGGG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CCAATATC

INDEX PRIMER	INDEX PRIMER SEQUENCE	EXPECTED INDEX SEQUENCE READ
P39-D3	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGCGTAGA</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TCTACGCA
P40-D4	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGCGTTAC</b> GTGACTGGAGTTCCAG ACGTGTGCTCTCCGATC-s-T-3'	GTAACCGA
P41-D5	5'-CAAGCAGAAGACGGCATAACGAGAT <b>ATGACGTC</b> GTGACTGGAGTTCCAG ACGTGTGCTCTCCGATC-s-T-3'	GACGTCAT
P42-D6	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GCTGTAA</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CTTACAGC
P43-D7	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGCAATGG</b> AGTACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TCCATTGC
P44-D8	5'-CAAGCAGAAGACGGCATAACGAGAT <b>ATCTCGCT</b> GTGACTGGAGTTCCAG ACGTGTGCTCTCCGATC-s-T-3'	AGCGAGAT
P45-D9	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GGCTATTG</b> TGACTGGAGTTCCAG ACGTGTGCTCTCCGATC-s-T-3'	CAATAGCC
P46-D10	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GGTGTCTT</b> GTGACTGGAGTTCCAG ACGTGTGCTCTCCGATC-s-T-3'	AAGACACC
P47-D11	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TCAACTGG</b> GTGACTGGAGTTCCAG GACGTGTGCTCTCCGATC-s-T-3'	CCAGTTGA
P48-D12	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CTTACCAC</b> GTGACTGGAGTTCCAG ACGTGTGCTCTCCGATC-s-T-3'	TGGTGAAG
P49-E1	5'-CAAGCAGAAGACGGCATAACGAGAT <b>ACGGCTT</b> GTGACTGGAGTTCCAG ACGTGTGCTCTCCGATC-s-T-3'	AAGACCGT
P50-E2	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TCTCGCAA</b> GTGACTGGAGTTCCAG ACGTGTGCTCTCCGATC-s-T-3'	TTGCGAGA
P51-E3	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GGAATTGC</b> TGACTGGAGTTCCAG GACGTGTGCTCTCCGATC-s-T-3'	GCAATTCC
P52-E4	5'-CAAGCAGAAGACGGCATAACGAGAT <b>ACGGATT</b> CTGACTGGAGTTCCAG ACGTGTGCTCTCCGATC-s-T-3'	GAATCCGT
P53-E5	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TTAAGCGG</b> GTGACTGGAGTTCCAG GACGTGTGCTCTCCGATC-s-T-3'	CCGCTTAA
P54-E6	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGCAGGTA</b> GTGACTGGAGTTCCAG GACGTGTGCTCTCCGATC-s-T-3'	TACCTGCA
P55-E7	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CAATCGAC</b> TGACTGGAGTTCCAG ACGTGTGCTCTCCGATC-s-T-3'	GTCGATTG
P56-E8	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GTGCCATA</b> GTGACTGGAGTTCCAG ACGTGTGCTCTCCGATC-s-T-3'	TATGGCAC
P57-E9	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGTTCCAG</b> GTGACTGGAGTTCCAG GACGTGTGCTCTCCGATC-s-T-3'	CTCGAACA
P58-E10	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGGAGTTG</b> TGACTGGAGTTCCAG GACGTGTGCTCTCCGATC-s-T-3'	CAACTCCA

## NEBNext Index/Universal Primer Mix (Cont.)

INDEX PRIMER	INDEX PRIMER SEQUENCE	EXPECTED INDEX SEQUENCE READ
P59-E11	5'-CAAGCAGAAGACGGCATAACGAGAT <b>ACGATGAC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GTCATCGT
P60-E12	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGATGTCCG</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GGACATCA
P61-F1	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGAACCTG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CAGGTTCA
P62-F2	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CTTCGTTG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GAACGAAG
P63-F3	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CTTCTGAG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CTCAGAAG
P64-F4	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGCTCATG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CATGAGCA
P65-F5	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AGTTCGTC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GACGAACT
P66-F6	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TAGCGTCT</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	AGACGCTA
P67-F7	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GGCGTTAT</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	ATAACGCC
P68-F8	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GGTGATTG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GAATCACC
P69-F9	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AACTTGCC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GGCAAGTT
P70-F10	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GCAAGATC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GATCTTGC
P71-F11	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TCCGATTG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CAATGCGA
P72-F12	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGTACACG</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GGTGACAA
P73-G1	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AGCTCCTA</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TAGGAGCT
P74-G2	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GCAATTCC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CGAATTGC
P75-G3	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CTTAGGAC</b> TGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GTCCTAAG
P76-G4	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GTCTTAAG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CTTAGGAC
P77-G5	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AACGTGGG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TCCACGTT
P78-G6	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CTGTGTTG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CAACACAG

INDEX PRIMER	INDEX PRIMER SEQUENCE	EXPECTED INDEX SEQUENCE READ
P79-G7	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GTAAAGGC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GCCTTAAC
P80-G8	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CACCTTAC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GTAAGGTG
P81-G9	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGGTAGCT</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	AGCTACCA
P82-G10	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CAGTGAAG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CTTCACTG
P83-G11	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GTTCAACCG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GGTTGAAC
P84-G12	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGGCTATC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GATAGCCA
P85-H1	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CTGGAGT</b> AGTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	TACTCCAG
P86-H2	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TCTCTCCG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GGAAGAGA
P87-H3	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TCTAACCG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GC GTTAGA
P88-H4	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GGTCAGAT</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	ATCTGACC
P89-H5	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CTCTGGTT</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	AACCAGAG
P90-H6	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGTGGTAC</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GTACCACA
P91-H7	5'-CAAGCAGAAGACGGCATAACGAGAT <b>CCTATACCG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	GGTATAGG
P92-H8	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TTCTCTCG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CGAGAGAA
P93-H9	5'-CAAGCAGAAGACGGCATAACGAGAT <b>GTATGCTG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CAGCATAC
P94-H10	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AAGTCGAG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CTCGACTT
P95-H11	5'-CAAGCAGAAGACGGCATAACGAGAT <b>AACCGAAG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CTTCGGTT
P96-H12	5'-CAAGCAGAAGACGGCATAACGAGAT <b>TGTTGTGG</b> GTGACTGGAGTTCA GACGTGTGCTCTCCGATC-s-T-3'	CCACAACA

## Quality Control Assays

**16-Hour Incubation:** 50  $\mu$ l reactions containing 1  $\mu$ l NEBNext Index/Universal Primer Mix and 1  $\mu$ g of HindIII digested Lambda DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50  $\mu$ l reactions containing NEBNext Index/Universal Primer Mix for Illumina and 1  $\mu$ g of T3 DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

**Endonuclease Activity:** Incubation of a 50  $\mu$ l reaction containing 1  $\mu$ l NEBNext Index/Universal Primer Mix with 1  $\mu$ g of  $\phi$ X174 RF I supercoiled DNA for 4 hours at 37°C results in less than 10% conversion to RF II (nicked molecules) as determined by agarose gel electrophoresis.

**RNase Activity:** Incubation of a 10  $\mu$ l reaction containing 1  $\mu$ l NEBNext Index/Universal Primer Mix with 40 ng of RNA transcript for 16 hours at 37°C resulted in no detectable degradation of RNA as determined by gel electrophoresis.

**Phosphatase Activity:** Incubation of NEBNext Index/Universal Primer Mix in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM  $MgCl_2$ ) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**Lot Controlled**

## Revision History:

REVISION #	DESCRIPTION	DATE
1.0	N/A	3/16
1.1	Corrected note on page 6 and 7 from "Do not use reuse primer if the seal has been previously pierced to avoid contamination with other indexed primers." to "Do not reuse primer if the seal has been previously pierced to avoid contamination with other indexed primers."	4/16
1.2	Updated Figure 1.1 page 7.	5/16
1.3	Clarified primer concentration. Corrected primer sequences 5' end.	6/16





## USA

New England Biolabs, Inc.  
240 County Road  
Ipswich, MA 01938-2723  
Telephone: (978) 927-5054  
Toll Free: (USA Orders) 1-800-632-5227  
Toll Free: (USA Tech) 1-800-632-7799  
Fax: (978) 921-1350  
e-mail: [info@neb.com](mailto:info@neb.com)  
[www.neb.com](http://www.neb.com)

## CANADA

New England Biolabs, Ltd.  
Telephone: (905) 665-4632  
Toll Free: 1-800-387-1095  
Fax: (905) 665-4635  
Fax Toll Free: 1-800-563-3789  
e-mail: [info.ca@neb.com](mailto:info.ca@neb.com)  
[www.neb.ca](http://www.neb.ca)

## CHINA, PEOPLE'S REPUBLIC

New England Biolabs (Beijing), Ltd.  
Telephone: 010-82378265/82378266  
Fax: 010-82378262  
e-mail: [info@neb-china.com](mailto:info@neb-china.com)  
[www.neb-china.com](http://www.neb-china.com)

## FRANCE

New England Biolabs France  
Free Call: 0800-100-632  
Free Fax: 0800-100-610  
e-mail: [info.fr@neb.com](mailto:info.fr@neb.com)  
[www.neb-online.fr](http://www.neb-online.fr)

## GERMANY & AUSTRIA

New England Biolabs GmbH  
Telephone: +49/(0)69/305 23140  
Free Call: 0800/246 5227 (Germany)  
Free Call: 00800/246 52277 (Austria)  
Fax: +49/(0)69/305 23149  
Free Fax: 0800/246 5229 (Germany)  
e-mail: [info.de@neb.com](mailto:info.de@neb.com)  
[www.neb-online.de](http://www.neb-online.de)

## JAPAN

New England Biolabs Japan, Inc.  
Telephone: +81 (0)3 5669 6191  
Fax: +81 (0)3 5669 6192  
e-mail: [info.jp@neb.com](mailto:info.jp@neb.com)  
[www.nebj.jp](http://www.nebj.jp)

## SINGAPORE

New England Biolabs Pte. Ltd.  
Telephone: +65 63859623  
Fax: +65 63859617  
e-mail: [sales.sg@neb.com](mailto:sales.sg@neb.com)  
[www.neb.sg](http://www.neb.sg)

## UNITED KINGDOM

New England Biolabs (UK) Ltd.  
Telephone: (01462) 420616  
Call Free: 0800 318486  
Fax: (01462) 421057  
Fax Free: 0800 435682  
e-mail: [info.uk@neb.com](mailto:info.uk@neb.com)  
[www.neb.uk.com](http://www.neb.uk.com)