

LIBRARY PREPARATION

NEBNext[®] Ultra[™] Directional RNA Library Prep Kit for Illumina[®]

Instruction Manual

NEB #E7420S/L
24/96 reactions

 NEW ENGLAND
BioLabs[®] Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.



This product is covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please contact NEB's Global Business Development team at gbd@neb.com.

AGENCOURT[®], AMPURE[®] and RNACLEAN[®] are registered trademarks of Beckman Coulter, Inc.

E-GEL[®] is a registered trademark of Life Technologies, Inc.

SIZESELECT[™] is a trademark of Life Technologies, Inc.

BIOANALYZER[®] is a registered trademark of Agilent Technologies, Inc.

ILLUMINA[®] is a registered trademark of Illumina, Inc.

MINELUTE[®] and QIAQUICK[®] are registered trademarks of Qiagen.

IGEPAL[®] is a registered trademark of Rhodia Operations.





Table of Contents:

| | |
|---|----|
| The Library Prep Kit Includes..... | 2 |
| Required Materials Not Included..... | 2 |
| Applications..... | 3 |
| Chapter 1 | |
| Protocol for use with NEBNext Poly(A) mRNA Magnetic Isolation Module..... | 5 |
| Chapter 2 | |
| Protocol for use with NEBNext rRNA Depletion Kit (Human/Mouse/Rat)..... | 19 |
| Chapter 3 | |
| Protocol for use with Purified mRNA or Ribosome Depleted RNA..... | 35 |
| Chapter 4 | |
| Appendix A..... | 47 |
| Chapter 5 | |
| Troubleshooting Guide..... | 51 |
| Chapter 6 | |
| Frequently Asked Questions (FAQs)..... | 55 |
| NEBNext First Strand Synthesis Reaction Buffer..... | 57 |
| Random Primers..... | 58 |
| ProtoScript II Reverse Transcriptase..... | 59 |
| Murine RNase Inhibitor..... | 60 |
| NEBNext Second Strand Synthesis Enzyme Mix..... | 61 |
| NEBNext Second Strand Synthesis Reaction Buffer with dUTP Mix..... | 62 |
| NEBNext End Prep Enzyme Mix..... | 63 |
| NEBNext End Repair Reaction Buffer..... | 64 |
| Blunt/TA Ligase Master Mix..... | 65 |
| Nuclease-free Water..... | 66 |
| NEBNext Q5 Hot Start HiFi PCR Master Mix..... | 67 |
| NEBNext USER™ Enzyme..... | 68 |
| Revision History..... | 69 |

The Library Prep Kit Includes:

The volumes provided are sufficient for preparation of up to 24 reactions (NEB #E7420S) and 96 reactions (NEB #E7420L). (All reagents should be stored at -20°C).

- (pink) NEBNext First Strand Synthesis Reaction Buffer (5X)
 - (pink) Random Primers
 - (pink) ProtoScript II Reverse Transcriptase
 - (pink) Murine RNase Inhibitor 
 - (orange) NEBNext Second Strand Synthesis Enzyme Mix
 - (orange) NEBNext Second Strand Synthesis Reaction Buffer with dUTP Mix (10X)
 - (green) NEBNext End Prep Enzyme Mix
 - (green) NEBNext End Repair Reaction Buffer (10X)
 - (red) Blunt/TA Ligase Master Mix
- Nuclease-free water
- (blue) NEBNext Q5 Hot Start HiFi PCR Master Mix
 - (blue) NEBNext USER Enzyme

Required Materials Not Included:

NEBNext Poly(A) mRNA Magnetic Isolation Module (NEB #E7490) or NEBNext rRNA Depletion Kit (NEB #E6310)

NEBNext Singleplex (NEB #E7350) or NEBNext Multiplex (NEB #E7335, #E7500, #E7710, #E7730, #E6609, #E7600) Oligos for Illumina or customer supplied oligos

Magnetic Rack (Alpaqua, cat #A001322 or equivalent)

80% Ethanol (freshly prepared)

0.1X TE, pH 8.0

10 mM Tris-HCl, pH 7.5–8.0

10 mM NaCl (optional)

Agencourt[®] AMPure[®] XP Beads (Beckman Coulter, Inc. #A63881)

Actinomycin D (Sigma# A1410, dissolved in dimethylsulfoxide [DMSO] to 5 µg/ul). See page 9, 23 or 37 for details.

Additional required materials not included for use with NEBNext rRNA Depletion Kit (NEB #E6310): Agencourt RNAClean[®] XP (Beckman Coulter, Inc., Cat #A63987)

Applications:

The NEBNext Ultra RNA Directional Library Prep Kit for Illumina contains enzymes and buffers that are ideally suited for cDNA library preparation for next-generation sequencing. Each of these components must pass rigorous quality control standards and are lot controlled, both individually and as a set of reagents.

Lot Control: The lots provided in the NEBNext Ultra RNA Directional Library Prep Kit for Illumina are managed separately and are qualified by additional functional validation. Individual reagents undergo standard enzyme activity and quality control assays, and also meet stringent criteria in the additional quality controls listed on each individual component page.

Functionally Validated: Each set of reagents is functionally validated together through construction and sequencing of a transcriptome library on an Illumina sequencing platform.

For larger volume requirements, customized and bulk packaging is available by purchasing through the OEM/Bulks department at NEB. Please contact OEM@neb.com for further information.

1

Protocol for use with NEBNext Poly(A) mRNA Magnetic Isolation Module (NEB #E7490)

NEBNext Ultra Directional
RNA Library Prep Kit for Illumina
Instruction Manual

Please refer to revision history for a summary of protocol updates

Symbols



This is a point where you can safely stop the protocol and store the samples prior to proceeding to the next step in the protocol.



This caution sign signifies a step in the protocol that has two paths leading to the same end point but is dependent on a user variable, like the type of RNA input.



Colored bullets indicate the cap color of the reagent to be added

The protocol has been optimized using high quality Universal Human Reference Total RNA. For PolyA mRNA selection, high quality RNA with RIN score > 7 (measured by Agilent Bioanalyzer®) is required.

RNA Sample Recommendations

The RNA sample should be free of salts (e.g., Mg²⁺, or guanidinium salts) or organics (e.g., phenol and ethanol).

Treat the RNA sample with DNase I to remove all traces of DNA. Remove DNase I after treatment.

AMPure XP Beads are required throughout the protocol. Allow beads to reach room temperature prior to use.

Starting Material: Total RNA (100 ng–1 µg) quantified by Bioanalyzer.

The protocol is optimized for approximately 200 bp RNA inserts. To generate libraries with longer RNA insert sizes, refer to Appendix A for recommended fragmentation times and size selection conditions.

1.1. Preparation of First Strand Reaction Buffer and Random Primer Mix

Prepare the First Strand Synthesis Reaction Buffer and Random Primer Mix (2X) as follows in a nuclease-free tube:

| | |
|--|-------|
| ● (pink) NEBNext First Strand Synthesis Reaction Buffer (5X) | 8 µl |
| ● (pink) NEBNext Random Primers | 2 µl |
| Nuclease-free water | 10 µl |
| <hr/> | |
| Total Volume | 20 µl |

Note: Keep the mix on ice.

1.2 mRNA Isolation, Fragmentation and Priming Starting with Total RNA

1. Dilute the total RNA with nuclease-free water to a final volume of 50 μ l in a nuclease-free 0.2 ml PCR tube and keep on ice.
2. Aliquot 20 μ l of NEBNext Oligo d(T)₂₅ beads into a nuclease-free 0.2 ml PCR tube.
3. Wash the beads by adding 100 μ l of RNA Binding Buffer (2X) to the beads. Pipette the entire volume up and down 6 times to mix thoroughly.
4. Place the tubes on the magnetic rack at room temperature for 2 minutes.
5. Remove and discard all of the supernatant from the tube. Take care not to disturb the beads.
6. Remove the tube from the magnetic rack.
7. Repeat steps 3–6.
8. Resuspend the beads in 50 μ l of RNA Binding Buffer (2X) and add the 50 μ l of total RNA sample from Step 1.
9. Place the tube on a thermal cycler and close the lid. Heat the sample at 65°C for 5 minutes and hold at 4°C to denature the RNA and facilitate binding of the poly-A mRNA to the beads.
10. Remove the tube from the thermal cycler when the temperature reaches 4°C.
11. Resuspend the beads. Pipette up and down slowly 6 times to mix thoroughly.
12. Place the tube on the bench and incubate at room temperature for 5 minutes to allow the mRNA to bind to the beads.
13. Resuspend the beads. Pipette up and down slowly 6 times to mix thoroughly.
14. Incubate for 5 minutes at room temperature to allow the RNA to bind to the beads.
15. Place the tube on the magnetic rack at room temperature for 2 minutes to separate the poly-A mRNA bound to the beads from the solution.
16. Remove and discard all of the supernatant. Take care not to disturb the beads.
17. Remove the tube from the magnetic rack.
18. Wash the beads by adding 200 μ l of Wash Buffer to the tube to remove unbound RNA. Pipette the entire volume up and down 6 times to mix thoroughly.

19. Place the tube on the magnetic rack at room temperature for 2 minutes.
20. Remove and discard all of the supernatant from the tube. Take care not to disturb the beads.
21. Remove the tube from the magnetic rack.
22. Repeat steps 18–21.
23. Add 50 μ l of Tris Buffer to each tube. Gently pipette the entire volume up and down 6 times to mix thoroughly.
24. Place the tube on the thermal cycler. Close the lid and heat the samples at 80°C for 2 minutes, then hold at 25°C to elute the Poly-A mRNA from the beads.
25. Remove the tube from the thermal cycler when the temperature reaches 25°C.
26. Add 50 μ l of RNA Binding Buffer (2X) to the sample to allow the mRNA to re-bind to the beads. Gently pipette the entire volume up and down 6 times to mix thoroughly.
27. Incubate the tube at room temperature for 5 minutes.
28. Resuspend the beads. Pipette up and down slowly 6 times to mix thoroughly.
29. Incubate for 5 minutes at room temperature to allow the RNA to bind to the beads.
30. Place the tube on the magnetic rack at room temperature for 2 minutes.
31. Remove and discard all of the supernatant from the tube. Take care not to disturb the beads.
32. Remove the tube from the magnetic rack.
33. Wash the beads by adding 200 μ l of Wash Buffer. Gently pipette the entire volume up and down 6 times to mix thoroughly.
34. Place the tube on the magnetic rack at room temperature for 2 minutes.
35. Remove and discard all of the supernatant from the tube. Take care not to disturb the beads.

Note: It is important to remove all of the supernatant to successfully fragment the mRNA in the subsequent steps. Spin down the tube. Place the tube on the magnetic rack and with a 10 μ l tip remove all of the wash buffer. Caution: Do not disturb beads that contain the mRNA.

36. Remove the tube from the magnetic rack.



Note: For RNA insert sizes > 200 nt, refer to Chapter 5 (Appendix A) for recommended fragmentation times in Step 37.

37. Elute mRNA from the beads by adding 15.5 μl of the First Strand Synthesis Reaction Buffer and Random Primer mix (2X) prepared in Section 1.1 incubating the sample at 94°C for 15 minutes. Immediately, place the tubes on the magnetic rack.
38. Collect the purified mRNA by transferring 13.5 μl of the supernatant to a clean nuclease-free PCR Tube.
39. Place the tube on ice and proceed directly to first strand cDNA synthesis.

1.3 First Strand cDNA Synthesis




Dilute Actinomycin D stock solution (5 $\mu\text{g}/\mu\text{l}$) to 0.1 $\mu\text{g}/\mu\text{l}$ in nuclease free water for immediate use.

Note: Dilute solutions of Actinomycin D are very sensitive to light. In solution, Actinomycin D tends to adsorb to plastic and glass. For these reasons, unused dilute solutions should be discarded and not stored for further use. However, frozen aliquots of a concentrated stock solution (5 $\mu\text{g}/\mu\text{l}$) in DMSO are expected to be stable for at least a month at -20°C.

1. To the fragmented and primed mRNA (13.5 μl from above Section 1.2, Step 38) add the following components and mix by gentle pipetting:

| | |
|--|-------------------|
| ● (pink) Murine RNase Inhibitor | 0.5 μl |
| Actinomycin D (0.1 $\mu\text{g}/\mu\text{l}$) | 5 μl |
| ● (pink) ProtoScript II Reverse Transcriptase | 1 μl |
| <hr/> | |
| Final volume | 20 μl |

Note: If you are following recommendations in Appendix A, for longer RNA fragments, increase the incubation at 42°C from 15 minutes to 50 minutes in Step 2.

2.  Incubate the sample in a preheated thermal cycler (with the heated lid set at 105°C) as follows:
 - 10 minutes at 25°C
 - 15 minutes at 42°C
 - 15 minutes at 70°C
 - Hold at 4°C
3. Immediately, perform second strand cDNA synthesis reaction.

1.4 Perform Second Strand cDNA Synthesis

1. Add the following reagents to the First Strand Synthesis reaction (20 μ l):

| | |
|--|------------|
| Nuclease-free water | 48 μ l |
| ● (orange) Second Strand Synthesis Reaction Buffer (10X) | 8 μ l |
| ● (orange) Second Strand Synthesis Enzyme Mix | 4 μ l |
| <hr/> | |
| Total volume | 80 μ l |
2. Mix thoroughly by gentle pipetting.
3. Incubate in a thermal cycler for 1 hour at 16°C, with heated lid set at \leq 40°C.

1.5 Purify the Double-stranded cDNA Using 1.8X Agencourt AMPure XP Beads

1. Vortex AMPure XP Beads to resuspend.
2. Add 144 μ l (1.8X) of resuspended AMPure XP Beads to the second strand synthesis reaction (~80 μ l). Mix well on a vortex mixer or by pipetting up and down at least 10 times.
3. Incubate for 5 minutes at room temperature.
4. Quickly spin the tube in a microcentrifuge to collect any sample on the sides of the tube. Place the tube on an appropriate magnetic rack to separate beads from supernatant. After the solution is clear (about 5 minutes), carefully remove and discard the supernatant. Be careful not to disturb the beads that contain DNA targets.
5. Add 200 μ l of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
6. Repeat Step 5 once for a total of 2 washing steps.
7. Air dry the beads for 5 minutes while the tube is on the magnetic rack with lid open.

Caution: Do not overdry the beads. This may result in lower recovery of DNA target.

8. Remove the tube from the magnet. Elute the DNA target from the beads into 60 μ l 0.1X TE Buffer or 10 mM Tris-HCl. Mix well on a vortex mixer or by pipetting up and down. Quickly spin the tube and incubate for 2 minutes at room temperature. Place the tube in the magnetic rack until the solution is clear.
9. Remove 55.5 μ l of the supernatant and transfer to a clean nuclease-free PCR tube.



Note: If you need to stop at this point in the protocol samples can be stored at -20°C.

1.6 Perform End Prep of cDNA Library

- Mix the following components in a sterile nuclease-free tube:

| | |
|---|--------------|
| Purified double stranded cDNA (Step 9, Section 1.5) | 55.5 μ l |
| ● (green) NEBNext End Repair Reaction Buffer (10X) | 6.5 μ l |
| ● (green) NEBNext End Prep Enzyme Mix | 3 μ l |
| <hr/> | |
| Total volume | 65 μ l |
- Incubate the sample in a thermal cycler (with the heated lid set at 75°C) as follows:
 - 30 minutes at 20°C
 - 30 minutes at 65°C
 - Hold at 4°C
- Proceed immediately to Adaptor Ligation.

1.7 Perform Adaptor Ligation

Dilute the ● (red) NEBNext Adaptor* for Illumina (15 μ M) to 1.5 μ M with a 10-fold dilution (1:9) with 10 mM Tris-HCl and 10 mM NaCl for immediate use.

- Add the following components directly to the end prep reaction mixture. **(Caution: do not pre-mix the components to prevent adaptor-dimer formation):**

| | |
|------------------------------------|--------------|
| End Prep Reaction | 65 μ l |
| ● (red) Blunt/TA Ligase Master Mix | 15 μ l |
| Diluted NEBNext Adaptor* | 1 μ l |
| Nuclease-free Water | 2.5 μ l |
| <hr/> | |
| Total volume | 83.5 μ l |

*The adaptor is provided in NEBNext Singleplex (NEB #E7350) or NEBNext Multiplex (NEB #E7335, #E7500, #E7710, #E7730, #E6609 or #E7600) Oligos for Illumina.

- Mix by pipetting followed by a quick spin to collect all liquid from the sides of the tube.
- Incubate 15 minutes at 20°C in a thermal cycler. Turn off the heated lid on the thermal cycler.



A precipitate can form upon thawing of the NEBNext Q5 Hot Start HiFi PCR Master Mix. To ensure optimal performance, place the master mix at room temperature while purifying the ligation reaction. Once thawed, gently mix by inverting the tube several times.

1.8 Purify the Ligation Reaction Using AMPure XP Beads



Note: If you are selecting for larger size fragments (> 200 nt) follow the size selection recommendations in Appendix A, Chapter 4.

1. Add nuclease-free water to the ligation reaction to bring the reaction volume to 100 μ l. It is important to ensure the final volume is 100 μ l prior to adding AMPure XP Beads.

Note: X refers to the original sample volume of 100 μ l from the above step.

2. Add 100 μ l (1.0X) resuspended AMPure XP Beads and mix well on a vortex mixer or by pipetting up and down at least 10 times.
3. Incubate for 5 minutes at room temperature.
4. Quickly spin the tube in a microcentrifuge and place the tube on an appropriate magnetic rack to separate beads from the supernatant. After the solution is clear (about 5 minutes), discard the supernatant that contain unwanted fragments (**Caution: do not discard the beads**).
5. Add 200 μ l of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
6. Repeat Step 5 once for a total of 2 washing steps.
7. Briefly spin the tube, and put the tube back in the magnetic rack.
8. Completely remove the residual ethanol, and air dry beads for 5 minutes while the tube is on the magnetic rack with the lid open.

Caution: Do not overdry the beads. This may result in lower recovery of DNA target.

9. Remove the tube from the rack. Elute DNA target from the beads with 52 μ l 0.1X TE or 10 mM Tris-HCl. Mix well on a vortex mixer or by pipetting up and down, and incubate for 2 minutes at room temperature. Put the tube in the magnetic rack until the solution is clear.
10. Transfer 50 μ l supernatant to a clean PCR tube. Discard beads.
11. To the 50 μ l supernatant, add 50 μ l (1.0X) of the resuspended AMPure XP Beads and mix well on a vortex mixer or by pipetting up and down at least 10 times.
12. Incubate for 5 minutes at room temperature.
13. Quickly spin the tube in a microcentrifuge and place the tube on an appropriate magnetic rack to separate beads from the supernatant. After the solution is clear (about 5 minutes), discard the supernatant that contains unwanted fragments (**Caution: do not discard the beads**).

14. Add 200 μ l of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
15. Repeat Step 14 once for a total of 2 washing steps.
16. Briefly spin the tube, and put the tube back in the magnetic rack.
17. Completely remove the residual ethanol, and air dry beads for 5 minutes while the tube is on the magnetic rack with the lid open.

Caution: Do not overdry the beads. This may result in lower recovery of DNA target.

18. Remove the tube from the rack. Elute DNA target from the beads with 19 μ l 0.1X TE or 10 mM Tris-HCl. Mix well on a vortex mixer or by pipetting up and down, and incubate for 2 minutes at room temperature. Put the tube in the magnetic rack until the solution is clear.
19. Without disturbing the bead pellet, transfer 17 μ l of the supernatant to a clean PCR tube and proceed to PCR enrichment.

1.9 PCR Enrichment of Adaptor Ligated DNA



Follow Section 1.9A if you are using the following oligos (10 μ M primer):

NEBNext Singleplex Oligos for Illumina (NEB #E7350)

NEBNext Multiplex Oligos for Illumina (Set 1, NEB #E7335)

NEBNext Multiplex Oligos for Illumina (Set 2, NEB #E7500)

NEBNext Multiplex Oligos for Illumina (Set 3, NEB #E7710)

NEBNext Multiplex Oligos for Illumina (Set 4, NEB #E7730)

NEBNext Multiplex Oligos for Illumina (Dual Index Primers, NEB #E7600)

Follow Section 1.9B if you are using NEBNext Multiplex Oligos for Illumina (96 Index Primers, NEB #E6609).

1.9A PCR Library Enrichment

- To the cDNA (17 μ l) from Step 19 Section 1.8 add the following components and mix by gentle pipetting:

| | |
|---|-------------|
| ● (blue) NEBNext USER Enzyme | 3 μ l |
| ● (blue) NEBNext Q5 Hot Start HiFi PCR Master Mix | 25 μ l |
| ● (blue) Index (X) Primer/i7 Primer*, ** | 2.5 μ l |
| ● (blue) Universal PCR Primer/i5 Primer*, *** | 2.5 μ l |
| Total volume | 50 μ l |

* The primers are provided in NEBNext Singleplex (NEB #E7350) or Multiplex (NEB #E7335, #E7500, #E7710, #E7730, #E7600) Oligos for Illumina. For use with Dual Index Primers (NEB #E7600), look at the NEB #E7600 manual for valid barcode combinations and tips for setting up PCR reactions.

** For use with NEBNext Multiplex Oligos (NEB #E7335, #E7710, #E7730 or #E7500) use only one Index Primer per PCR reaction. For use with Dual Index Primers (NEB #E7600) use only one i7 Primer per reaction.

*** For use with Dual Index Primers (NEB #E7600) use only one i5 Primer per reaction.

- PCR Cycling Conditions

| CYCLE STEP | TEMP | TIME | CYCLES |
|----------------------|------|------------|------------|
| USER Digestion | 37°C | 15 minutes | 1 |
| Initial Denaturation | 98°C | 30 seconds | 1 |
| Denaturation | 98°C | 10 seconds | 12–15*, ** |
| Annealing/Extension | 65°C | 75 seconds | |
| Final Extension | 65°C | 5 minutes | 1 |
| Hold | 4°C | ∞ | |

* The number of PCR cycles should be adjusted based on RNA input. If 100 ng total RNA is the starting input, it is recommended to perform 15 cycles of PCR.

** It is important to limit the number of PCR cycles to avoid overamplification. If overamplification occurs, larger molecular weight products (> 500 bp) will appear on the Bioanalyzer trace.

- Proceed to Section 1.10 (Purify the PCR Reaction using Agencourt AMPure XP Beads).

1.9B PCR Library Enrichment

- To the cDNA (17 μ l) from Step 19 Section 1.8 add the following components and mix by gentle pipetting:

| | |
|---|------------|
| ● (blue) NEBNext USER Enzyme | 3 μ l |
| ● (blue) Index/ Universal Primer Mix* | 5 μ l |
| ● (blue) NEBNext Q5 Hot Start HiFi PCR Master Mix | 25 μ l |

| | |
|--------------|------------|
| Total volume | 50 μ l |
|--------------|------------|

* The primers are provided in NEBNext Multiplex Oligos for Illumina, NEB #E6609. Refer to NEB #E6609 manual for valid barcode combinations and tips for setting up PCR reactions.

- PCR Cycling Conditions

| CYCLE STEP | TEMP | TIME | CYCLES |
|----------------------|------|------------|------------|
| USER Digestion | 37°C | 15 minutes | 1 |
| Initial Denaturation | 98°C | 30 seconds | 1 |
| Denaturation | 98°C | 10 seconds | 12–15*, ** |
| Annealing/Extension | 65°C | 75 seconds | |
| Final Extension | 65°C | 5 minutes | 1 |
| Hold | 4°C | ∞ | |

* The number of PCR cycles should be adjusted based on RNA input. If 100 ng total RNA is the starting input, it is recommended to perform 15 cycles of PCR.

** It is important to limit the number of PCR cycles to avoid overamplification. If overamplification occurs, larger molecular weight products (> 500 bp) will appear on the Bioanalyzer trace.

- Proceed to Section 1.10 (Purify the PCR Reaction using Agencourt AMPure XP Beads).

1.10 Purify the PCR Reaction using Agencourt AMPure XP Beads

Note: X refers to the original sample volume from the above step.

1. Vortex Agencourt AMPure XP Beads to resuspend.
2. Add 45 μl (0.9X) of resuspended Agencourt AMPure XP Beads to the PCR reaction (~ 50 μl). Mix well on a vortex mixer or by pipetting up and down at least 10 times.
3. Incubate for 5 minutes at room temperature.
4. Quickly spin the tube in a microcentrifuge and place the tube on an appropriate magnetic rack to separate beads from the supernatant. After the solution is clear (about 5 minutes), carefully remove and discard the supernatant. Be careful not to disturb the beads that contain DNA targets.
5. Add 200 μl of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
6. Repeat Step 5 once for a total of 2 washing steps.
7. Air dry the beads for 5 minutes while the tube is on the magnetic rack with the lid open.

Caution: Do not overdry the beads. This may result in lower recovery of DNA target.


8. Remove the tube from the rack. Elute the DNA target from the beads into 23 μl 0.1X TE. Mix well on a vortex mixer or by pipetting up and down, quickly spin the tube in a microcentrifuge and incubate for 2 minutes at room temperature. Place it in the magnetic rack until the solution is clear.
9. Transfer 20 μl of the supernatant to a clean PCR tube, and store at -20°C .

1.11 Assess library quality on a Bioanalyzer® (Agilent High Sensitivity Chip)

1. Dilute 2–3 μl of the library in 10 mM Tris or 0.1X TE.
2. Run 1 μl in a DNA High Sensitivity Chip
3. Check that the electropherogram shows a narrow distribution with a peak size approximately 300 bp.

Note: If a peak at ~ 80 bp (primers) or 128 bp (adaptor-dimer) is shown in the Bioanalyzer traces; Bring up the sample volume (Step 9, Section 1.10) to 50 μl exactly with nuclease-free water and repeat the AMPure XP Bead clean up step (Section 1.10).

Figure 1.1: Example of RNA library size distribution on a Bioanalyzer.



2

Protocol for use with NEBNext rRNA Depletion Kit (Human/Mouse/Rat) (NEB #E6310)

NEBNext Ultra Directional
RNA Library Prep Kit for Illumina
Instruction Manual

Note: *There is a formulation change to one of the components of the kit (NEB #E6318: NEBNext RNase H). The protocol and use of this component have not changed.*

Symbols



This is a point where you can safely stop the protocol and store the samples prior to proceeding to the next step in the protocol.



This caution sign signifies a step in the protocol that has two paths leading to the same end point but is dependent on a user variable, like the type of RNA input.



Colored bullets indicate the cap color of the reagent to be added

This protocol has been optimized using Universal Human Reference Total RNA.

RNA Sample Recommendations

The RNA sample should be free of salts (e.g., Mg^{2+} , or guanidinium salts), divalent cation chelating agents (e.g. EDTA, EGTA, citrate), or organics (e.g., phenol and ethanol).

Treat the RNA sample with DNase I to remove all traces of DNA. Remove DNaseI after treatment.

Typical Yield of rRNA-depleted RNA from a Reaction

The actual yield is dependent on the quality of the input RNA, the rRNA content of the sample, and the method used to purify the rRNA-depleted RNA. Recoveries of 3%–10% of the input RNA are typical.

RNA Input

100 ng to 1 μ g total RNA in up to 12 μ l total volume.

Note: The NEBNext rRNA Depletion Kit can be used with as low as 10 ng total RNA, however, for RNAseq samples we recommend using total RNA inputs 100 ng–1 μ g to increase library complexity and reduce sequencing duplication rates.



Assess quality of the input RNA by running input RNA on an Agilent Bioanalyzer to determine the RIN number (RIN). Highly degraded samples (RIN #1–2) (FFPE) or partially degraded samples (RIN #2–7) will require different fragmentation times (Section 2.5).

2.1 Hybridize the Probes to the RNA

1. Prepare a RNA/Probe master mix as follows:

| | |
|---------------------------------|-----------|
| NEBNext rRNA Depletion Solution | 1 μ l |
| Probe Hybridization Buffer | 2 μ l |
| <hr/> | |
| Total Volume | 3 μ l |

2. Add 3 μ l of the above mix to 12 μ l total RNA sample.
3. Mix by pipetting up and down.
4. Spin down briefly in a tabletop centrifuge, and immediately proceed to the next step.
5. Place samples in a thermocycler, and run the following program with the heated lid set at 105°C. This will take approximately 15–20 minutes to complete:
 - 2 minutes at 95°C
 - 0.1°C/sec at 95–22°C
 - 5 minutes hold at 22°C
6. Spin down the samples in a tabletop centrifuge, and place on ice. Proceed immediately to the next step.

2.2 RNase H Digestion

1. On ice, prepare a master mix according to the following, and mix by pipetting up and down; use immediately.

| | |
|-------------------------|-----------|
| NEBNext RNase H | 2 μ l |
| RNase H Reaction Buffer | 2 μ l |
| Nuclease-free Water | 1 μ l |
| <hr/> | |
| Total Volume | 5 μ l |

2. Add 5 μ l of the above mix to the RNA sample from Step 6 in Section 2.1.
3. Mix by pipetting up and down.
4. Place samples in a thermocycler (with lid at 40°C) and incubate at 37°C for 30 minutes.
5. Spin down the samples in a tabletop centrifuge, and place on ice. Proceed immediately to the next step to prevent non-specific degradation of RNA.

2.3 DNase I Digestion

1. On ice, prepare a DNase I Digestion Master Mix according to the following table, and mix by pipetting up and down; use immediately

| | |
|-------------------------|--------------|
| DNase I Reaction Buffer | 5 μ l |
| DNase I (RNase-free) | 2.5 μ l |
| Nuclease-free Water | 22.5 μ l |
| <hr/> | |
| Total volume | 30 μ l |

2. Add 30 μ l of the above mix to the RNA sample from Step 5 in Section 2.2, and mix by pipetting up and down.
3. Place samples in a thermocycler (with lid at 40°C) and incubate at 37°C for 30 minutes.
4. Spin down the samples in a tabletop centrifuge, and place on ice. Proceed immediately to the next step.

2.4 RNA Purification after rRNA Depletion Using Agencourt RNAClean XP

1. Add 2.2X (110 μ l) Agencourt RNAClean XP Beads to the RNA sample from the previous section (Section 2.3, Step 4) and mix by pipetting up and down.
2. Incubate samples on ice for 15 minutes.
3. Place the tube on an appropriate magnetic rack to separate beads from the supernatant.
4. When the solution is clear (about 5 minutes), discard the supernatant.
5. Add 200 μ l of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
6. Repeat Step 5 once for a total of 2 washes.
7. Briefly spin the tube, and put the tube back in the magnetic rack.
8. Completely remove the residual ethanol, and air dry beads.
9. Remove the tube from the magnetic rack. Elute RNA from the beads with 8 μ l nuclease-free water.
10. Mix well by pipetting up and down, and put the tube in the magnetic rack until the solution is clear.
11. Transfer 6 μ l of the supernatant to a clean PCR tube.
12. Place the sample on ice and proceed to Section 2.5.

2.5 RNA Fragmentation, Priming and First Strand cDNA Synthesis



RNA fragmentation is only required for intact or partially degraded RNA. Recommended fragmentation times can be found in Table 2.1. Follow protocol in 2.5A to set up the reaction. For highly degraded RNA (FFPE Samples) which do not require fragmentation proceed to Step 2.5B.

2.5A RNA Fragmentation and Priming Starting from Intact or Partially Degraded RNA:

- Set up the following reaction and mix by gentle pipetting:

| | |
|--|------------|
| Ribosome depleted RNA | 5 μ l |
| ● (pink) NEBNext First Strand Synthesis Reaction Buffer (5X) | 4 μ l |
| ● (pink) Random Primers | 1 μ l |
| Final volume | 10 μ l |

-  Incubate the sample at 94°C following the recommendations in Table 2.1 for fragments sizes ~200 nt.

Table 2.1. Suggested fragmentation times based on RIN number of RNA input.

| RNA TYPE | RIN | FRAG. TIME |
|------------------------|-----|------------------|
| Intact RNA | >7 | 15 min. at 94°C |
| Partially Degraded RNA | 2–6 | 7–8 min. at 94°C |

Refer to Appendix A for fragmentation conditions if you are preparing libraries with large inserts (> 200 bp). Conditions in Appendix A only apply for intact RNA.

- Transfer the tube to ice.

First Strand cDNA Synthesis




Dilute Actinomycin D stock solution (5 μ g/ μ l) to 0.1 μ g/ μ l in nuclease free water for immediate use.

Note: Dilute solutions of Actinomycin D are very sensitive to light. In solution, Actinomycin D tends to adsorb to plastic and glass. For these reasons, unused dilute solutions should be discarded and not stored for further use. However, frozen aliquots of a concentrated stock solution (5 μ g/ μ l) in DMSO are expected to be stable for at least a month at -20°C.

4. To the fragmented and primed mRNA from Step 3 in Section 2.5A (10 μ l), add the following components and mix by gentle pipetting:

| | |
|---|-------------|
| ● (pink) Murine RNase Inhibitor | 0.5 μ l |
| Actinomycin D (0.1 μ g/ μ l) | 5 μ l |
| ● (pink) ProtoScript II Reverse Transcriptase | 1 μ l |
| Nuclease free water | 3.5 μ l |
| <hr/> | |
| Final volume | 20 μ l |

Note: If you are following recommendations in Appendix A, for longer RNA fragments, increase the incubation at 42°C from 15 minutes to 50 minutes in Step 5.

5.  Incubate the sample in a preheated thermal cycler (with the heated lid set at 105°C) as follows:
- 10 minutes at 25°C
 - 15 minutes at 42°C
 - 15 minutes at 70°C
 - Hold at 4°C
6. Proceed directly to Second Strand cDNA Synthesis, Section 2.6.

2.5B Priming of Highly Degraded RNA (FFPE) which has a RIN ≤ 2 and does not Require Fragmentation:

1. Set up the following priming reaction and mix by gentle pipetting:
- | | |
|-------------------------|-----------|
| Ribosome depleted RNA | 5 μ l |
| ● (pink) Random Primers | 1 μ l |
| <hr/> | |
| Final volume | 6 μ l |
2. Incubate the sample in a preheated thermal cycler as follows:
5 minutes at 65°C, with heated lid set at 105°C. Hold at 4°C.
3. Transfer the tube directly to ice.

First Strand cDNA Synthesis




Dilute Actinomycin D stock solution (5 μ g/ μ l) to 0.1 μ g/ μ l in nuclease free water for immediate use.

Note: Dilute solutions of Actinomycin D are very sensitive to light. In solution, Actinomycin D tends to adsorb to plastic and glass. For these reasons, unused dilute solutions should be discarded and not stored for further use. However, frozen aliquots of a concentrated stock solution (5 μ g/ μ l) in DMSO are expected to be stable for at least a month at -20°C.

4. To the primed RNA from Step 3 Section 2.5B (6 μ l) add the following components and mix by gentle pipetting:

| | |
|--|-------------|
| ● (pink) NEBNext First Strand Synthesis Reaction Buffer (5X) | 4 μ l |
| ● (pink) Murine RNase Inhibitor | 0.5 μ l |
| Actinomycin D (0.1 μ g/ μ l) | 5 μ l |
| ● (pink) ProtoScript II Reverse Transcriptase | 1 μ l |
| Nuclease free water | 3.5 μ l |
| <hr/> | |
| Final volume | 20 μ l |

Note: If you are following recommendations in Appendix A, for longer RNA fragments, increase the incubation at 42°C from 15 minutes to 50 minutes in Step 5.

5.  Incubate the sample in a preheated thermal cycler (with the heated lid set at 105°C) as follows:
- 10 minutes at 25°C
 - 15 minutes at 42°C
 - 15 minutes at 70°C
 - Hold at 4°C
6. Proceed directly to Second Strand cDNA Synthesis, Section 2.6.

2.6 Perform Second Strand cDNA Synthesis

1. Add the following reagents to the First Strand Synthesis reaction (20 μ l):
- | | |
|--|------------|
| Nuclease-free water | 48 μ l |
| ● (orange) Second Strand Synthesis Reaction Buffer (10X) | 8 μ l |
| ● (orange) Second Strand Synthesis Enzyme Mix | 4 μ l |
| <hr/> | |
| Total volume | 80 μ l |
2. Mix thoroughly by gentle pipetting.
3. Incubate in a thermal cycler for 1 hour at 16°C, with heated lid set at \leq 40°C.

2.7 Purify the Double-stranded cDNA Using 1.8X Agencourt AMPure XP Beads

- Vortex AMPure XP Beads to resuspend.
- Add 144 μ l (1.8X) of resuspended AMPure XP Beads to the second strand synthesis reaction (~80 μ l). Mix well on a vortex mixer or by pipetting up and down at least 10 times.

3. Incubate for 5 minutes at room temperature.
4. Quickly spin the tube in a microcentrifuge to collect any sample on the sides of the tube. Place the tube on an appropriate magnetic rack to separate beads from supernatant. After the solution is clear (about 5 minutes), carefully remove and discard the supernatant. Be careful not to disturb the beads that contain DNA targets.
5. Add 200 μl of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
6. Repeat Step 5 once for a total of 2 washing steps.
7. Air dry the beads for 5 minutes while the tube is on the magnetic rack with lid open.

Caution: Do not over dry the beads. This may result in lower recovery of DNA target.

8. Remove the tube from the magnet. Elute the DNA target from the beads into 60 μl 0.1X TE Buffer or 10 mM Tris-HCl. Mix well on a vortex mixer or by pipetting up and down. Quickly spin the tube and incubate for 2 minutes at room temperature. Place the tube in the magnetic rack until the solution is clear.
9. Remove 55.5 μl of the supernatant and transfer to a clean nuclease free PCR tube.



Note: If you need to stop at this point in the protocol samples can be stored at -20°C .

2.8 Perform End Prep of cDNA Library

1. Mix the following components in a sterile nuclease free tube:

| | |
|---|--------------------|
| Purified double-stranded cDNA (Step 9, Section 2.7) | 55.5 μl |
| ● (green) NEBNext End Repair Reaction Buffer (10X) | 6.5 μl |
| ● (green) NEBNext End Prep Enzyme Mix | 3 μl |
| <hr/> | |
| Total volume | 65 μl |
2. Incubate the sample in a thermal cycler (with the heated lid set at 75°C) as follows:
30 minutes at 20°C
30 minutes at 65°C
Hold at 4°C
3. Proceed immediately to Adaptor Ligation.

2.9 Perform Adaptor Ligation



Dilute the ● (red) NEBNext adaptor* prior to setting up the ligation reaction.

| INPUT RNA | DILUTION REQUIRED |
|---------------------|--|
| 100 ng | 30 fold dilution in 10 mM Tris-HCl with 10 mM NaCl |
| > 100 ng up to 1 µg | 10 fold dilution in 10 mM Tris-HCl with 10 mM NaCl |

1. Add the following components directly to the End Prep Reaction (**Caution: Do not pre-mix the components to prevent adaptor-dimer formation**):

| | |
|------------------------------------|---------|
| End Prep Reaction | 65 µl |
| ● (red) Blunt/TA Ligase Master Mix | 15 µl |
| Diluted NEBNext Adaptor* | 1 µl |
| Nuclease-free Water | 2.5 µl |
| <hr/> | |
| Total volume | 83.5 µl |

*The adaptor is provided in NEBNext Singleplex (NEB #E7350) or NEBNext Multiplex (NEB #E7335, #E7500, #E7710, #E7730, #E6609 or #E7600) Oligos for Illumina.

2. Mix by pipetting followed by a quick spin to collect all liquid from the sides of the tube.
3. Incubate 15 minutes at 20°C in a thermal cycler.



A precipitate can form upon thawing of the NEBNext Hot Start HiFi PCR Master Mix. To ensure optimal performance, place the master mix at room temperature while purifying the ligation reaction. Once thawed, gently mix by inverting the tube several times.

2.10 Purify the Ligation Reaction Using AMPure XP Beads



Note: If you are selecting for larger size fragments (> 200 nt) follow the size selection recommendations in Appendix A, Chapter 4.

1. Add nuclease-free water to the ligation reaction to bring the reaction volume to 100 µl. It is important to ensure the final volume is 100 µl prior to adding AMPure XP Beads.

Note: X refers to the original sample volume of 100 µl from the above step.

2. Add 100 µl (1.0X) resuspended AMPure XP Beads and mix well on a vortex mixer or by pipetting up and down at least 10 times.

3. Incubate for 5 minutes at room temperature.
4. Quickly spin the tube in a microcentrifuge and place the tube on an appropriate magnetic rack to separate beads from the supernatant. After the solution is clear (about 5 minutes), discard the supernatant that contain unwanted fragments (**Caution: do not discard the beads**).
5. Add 200 μ l of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
6. Repeat Step 5 once for a total of 2 washing steps.
7. Briefly spin the tube, and put the tube back in the magnetic rack.
8. Completely remove the residual ethanol, and air dry beads for 5 minutes while the tube is on the magnetic rack with the lid open.
Caution: Do not over dry the beads. This may result in lower recovery of DNA target.
9. Remove the tube from the magnet. Elute DNA target from the beads with 52 μ l 0.1X TE or 10 mM Tris-HCl. Mix well on a vortex mixer or by pipetting up and down, incubate for 2 minutes at room temperature. Put the tube in the magnetic rack until the solution is clear.
10. Transfer the 50 μ l supernatant to a clean PCR tube. Discard beads.
11. To the 50 μ l supernatant, add 50 μ l (1.0X) of the resuspended AMPure XP Beads and mix well on a vortex mixer or by pipetting up and down at least 10 times.
12. Incubate for 5 minutes at room temperature.
13. Quickly spin the tube in a microcentrifuge and place the tube on an appropriate magnetic rack to separate beads from the supernatant. After the solution is clear (about 5 minutes), discard the supernatant that contains unwanted fragments (**Caution: do not discard the beads**).
14. Add 200 μ l of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
15. Repeat Step 14 once for a total of 2 washing steps.
16. Briefly spin the tube, and put the tube back in the magnetic rack.
17. Completely remove the residual ethanol, and air dry beads for 5 minutes while the tube is on the magnetic rack with the lid open.

Caution: Do not over dry the beads. This may result in lower recovery of DNA target.

18. Remove the tube from the magnet. Elute DNA target from the beads with 19 μ l 0.1X TE or 10 mM Tris-HCl. Mix well on a vortex mixer or by pipetting up and down, incubate for 2 minutes at room temperature. Put the tube in the magnetic rack until the solution is clear.
19. Without disturbing the bead pellet, transfer 17 μ l of the supernatant to a clean PCR tube and proceed to PCR enrichment.

2.11 PCR Enrichment of Adaptor Ligated DNA



Follow Section 2.11A if you are using the following oligos (10 μ M primer):

NEBNext Singleplex Oligos for Illumina (NEB #E7350)

NEBNext Multiplex Oligos for Illumina (Set 1, NEB #E7335)

NEBNext Multiplex Oligos for Illumina (Set 2, NEB #E7500)

NEBNext Multiplex Oligos for Illumina (Set 3, NEB #E7710)

NEBNext Multiplex Oligos for Illumina (Set 4, NEB #E7730)

NEBNext Multiplex Oligos for Illumina (Dual Index Primers, NEB #E7600)

Follow Section 2.11B if you are using NEBNext Multiplex Oligos for Illumina (96 Index Primers, NEB #E6609).

2.11A PCR Library Enrichment

1. To the cDNA (17 μ l) from Step 19 Section 2.10 add the following components and mix by gentle pipetting:

| | |
|---|-------------|
| ● (blue) NEBNext USER Enzyme | 3 μ l |
| ● (blue) NEBNext Q5 Hot Start HiFi PCR Master Mix | 25 μ l |
| ● (blue) Index (X) Primer/i7 Primer*, ** | 2.5 μ l |
| ● (blue) Universal PCR Primer/i5 Primer*, *** | 2.5 μ l |
| <hr/> | |
| Total volume | 50 μ l |

* The primers are provided in NEBNext Singleplex (NEB #E7350) or Multiplex (NEB #E7335, #E7500, #E7710, #E7730, #E7600) Oligos for Illumina. For use with Dual Index Primers (NEB #E7600), look at the NEB #E7600 manual for valid barcode combinations and tips for setting up PCR reactions.

** For use with NEBNext Multiplex Oligos (NEB #E7335, #E7500, #E7710 or #E7730) use only one Index Primer per PCR reaction. For use with Dual Index Primers (NEB #E7600) use only one i7 Primer per reaction.

*** For use with Dual Index Primers (NEB #E7600) use only one i5 Primer per reaction.

2. PCR Cycling Conditions

| CYCLE STEP | TEMP | TIME | CYCLES |
|----------------------|------|------------|------------|
| USER Digestion | 37°C | 15 minutes | 1 |
| Initial Denaturation | 98°C | 30 seconds | 1 |
| Denaturation | 98°C | 10 seconds | 12–15*, ** |
| Annealing/Extension | 65°C | 75 seconds | |
| Final Extension | 65°C | 5 minutes | 1 |
| Hold | 4°C | ∞ | |

* The number of PCR cycles should be adjusted based on RNA input. If 100 ng total RNA is the starting input, it is recommended to perform 15 cycles of PCR. However, optimization of PCR cycle number may be required to avoid overamplification.

** It is important to limit the number of PCR cycles to avoid overamplification. If overamplification occurs, larger molecular weight products (> 500 bp) will appear on the Bioanalyzer trace.

3. Proceed to Section 2.12 (Purify the PCR Reaction using Agencourt AMPure XP Beads).

2.11B PCR Library Enrichment

- To the cDNA (17 μ l) from Step 19 Section 2.10 add the following components and mix by gentle pipetting:

| | |
|---|-----------------------------|
| ● (blue) NEBNext USER Enzyme | 3 μ l |
| ● (blue) NEBNext Q5 Hot Start HiFi PCR Master Mix | 25 μ l |
| ● (blue) Index/Universal Primer Mix* | 5 μ l |
| Total volume | 50 μl |

* The primers are provided in NEBNext Multiplex Oligos for Illumina, NEB #E6609. Refer to NEB #E6609 manual for valid barcode combinations and tips for setting up PCR reactions.

- PCR Cycling Conditions

| CYCLE STEP | TEMP | TIME | CYCLES |
|----------------------|------|------------|------------|
| USER Digestion | 37°C | 15 minutes | 1 |
| Initial Denaturation | 98°C | 30 seconds | 1 |
| Denaturation | 98°C | 10 seconds | 12–15*, ** |
| Annealing/Extension | 65°C | 75 seconds | |
| Final Extension | 65°C | 5 minutes | 1 |
| Hold | 4°C | ∞ | |

* The number of PCR cycles should be adjusted based on RNA input. If 100 ng total RNA is the starting input, it is recommended to perform 15 cycles of PCR.

** It is important to limit the number of PCR cycles to avoid overamplification.

If overamplification occurs, larger molecular weight products (> 500 bp) will appear on the Bioanalyzer trace.

- Proceed to Section 2.12 (Purify the PCR Reaction using Agencourt AMPure XP Beads).

2.12 Purify the PCR Reaction using Agencourt AMPure XP Beads

Note: X refers to the original sample volume from the above step.

1. Vortex Agencourt AMPure XP Beads to resuspend.
2. Add 45 μl (0.9X) of resuspended Agencourt AMPure XP Beads to the PCR reaction ($\sim 50 \mu\text{l}$). Mix well on a vortex mixer or by pipetting up and down at least 10 times.
3. Incubate for 5 minutes at room temperature.
4. Quickly spin the tube in a microcentrifuge and place the tube on an appropriate magnetic rack to separate beads from the supernatant. After the solution is clear (about 5 minutes), carefully remove and discard the supernatant. Be careful not to disturb the beads that contain DNA targets.
5. Add 200 μl of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
6. Repeat Step 5 once for a total of 2 washing steps.
7. Air dry the beads for 5 minutes while the tube is on the magnetic rack with the lid open.

Caution: Do not over dry the beads. This may result in lower recovery of DNA target.


8. Remove the tube from the magnetic rack. Elute the DNA target from the beads into 23 μl 0.1X TE. Mix well on a vortex mixer or by pipetting up and down, quickly spin the tube in a microcentrifuge and incubate for 2 minutes at room temperature. Place the tube in the magnetic rack until the solution is clear.
9. Transfer 20 μl of the supernatant to a clean PCR tube, and store at -20°C .

2.13 Assess library quality on a Bioanalyzer (Agilent High Sensitivity Chip).

1. Dilute 2–3 μl of the library in 10 mM Tris or 0.1X TE.
2. Run 1 μl in a DNA High Sensitivity chip
3. Check that the electropherogram shows a narrow distribution with a peak size approximately 300 bp.

Note: If a peak at ~ 80 bp (primers) or 128 bp (adaptor-dimer) is shown in the Bioanalyzer traces; Bring up the sample volume (Step 9, Section 2.12) to 50 μl exactly with nuclease-free water and repeat the AMPure XP Bead clean up step (Section 2.12).

Figure 2.1: Example of RNA library size distribution on a Bioanalyzer.





3

Protocol for use with Purified mRNA or Ribosome Depleted RNA

NEBNext Ultra Directional
RNA Library Prep Kit for Illumina
Instruction Manual

Please refer to revision history for a summary of protocol updates

Symbols



This is a point where you can safely stop the protocol and store the samples prior to proceeding to the next step in the protocol.



This caution sign signifies a step in the protocol that has two paths leading to the same end point but is dependent on a user variable, like the type of RNA input.



Colored bullets indicate the cap color of the reagent to be added

This protocol has been optimized using Universal Human Reference Total RNA.

3.1 RNA Fragmentation, Priming and First Strand cDNA Synthesis



RNA fragmentation is only required for intact or partially degraded RNA. Recommended fragmentation times can be found in Table 3.1. Follow protocol in 3.1A to set up the reaction. For highly degraded RNA (FFPE Samples) which do not require fragmentation proceed to Step 3.1B.

3.1A RNA Fragmentation and Priming Starting from Intact or Partially Degraded RNA:

- Set up the following reaction and mix by gentle pipetting:

| | |
|--|------------|
| Purified mRNA or Ribosome depleted RNA (10-100 ng) | 5 μ l |
| ● (pink) NEBNext First Strand Synthesis Reaction Buffer (5X) | 4 μ l |
| ● (pink) Random Primers | 1 μ l |
| Final volume | 10 μ l |

- Incubate the sample at 94°C following the recommendations in Table 3.1 for fragments sizes ~200 nt.

Table 3.1. Suggested fragmentation times based on RIN number of RNA input.

| RNA TYPE | RIN | FRAG. TIME |
|------------------------|-----|------------------|
| Intact RNA | >7 | 15 min. at 94°C |
| Partially Degraded RNA | 2–6 | 7–8 min. at 94°C |

Refer to Appendix A for fragmentation conditions if you are preparing libraries with large inserts (> 200 bp). Conditions in Appendix A only apply for intact RNA.

- Transfer the tube to ice.

First Strand cDNA Synthesis




Dilute Actinomycin D stock solution (5 µg/µl) to 0.1 µg/µl in nuclease free water for immediate use.

Note: Dilute solutions of Actinomycin D are very sensitive to light. In solution, Actinomycin D tends to adsorb to plastic and glass. For these reasons, unused dilute solutions should be discarded and not stored for further use. However, frozen aliquots of a concentrated stock solution (5 µg/µl) in DMSO are expected to be stable for at least a month at -20°C.

4. To the fragmented and primed mRNA from Step 3 in Section 3.1A (10 µl) add the following components and mix by gentle pipetting:

| | |
|---|--------|
| ● (pink) Murine RNase Inhibitor | 0.5 µl |
| Actinomycin D (0.1 µg/µl) | 5 µl |
| ● (pink) ProtoScript II Reverse Transcriptase | 1 µl |
| Nuclease free water | 3.5 µl |
| <hr/> | |
| Final volume | 20 µl |

Note: If you are following recommendations in Appendix A, for longer RNA fragments, increase the incubation at 42°C from 15 minutes to 50 minutes in Step 5.

5.  Incubate the sample in a preheated thermal cycler (with the heated lid set at 105°C) as follows:

10 minutes at 25°C
 15 minutes at 42°C
 15 minutes at 70°C
 Hold at 4°C

6. Proceed directly to Second Strand cDNA Synthesis, Section 3.2.

3.1B Priming of Highly Degraded RNA (FFPE) which has a RIN ≤ 2 and does not Require Fragmentation:

1. Set up the following priming reaction and mix by gentle pipetting:

| | |
|--|------|
| Purified mRNA or Ribosome depleted RNA (10-100 ng) | 5 µl |
| ● (pink) Random Primers | 1 µl |
| <hr/> | |
| Final volume | 6 µl |

2. Incubate the sample in a preheated thermal cycler as follows: 5 minutes at 65°C, with heated lid set at 105°C. Hold at 4°C.
3. Transfer the tube directly to ice.

First Strand cDNA Synthesis


-  Dilute Actinomycin D stock solution (5 µg/µl) to 0.1 µg/µl in nuclease free water for immediate use.

Note: Dilute solutions of Actinomycin D are very sensitive to light. In solution, Actinomycin D tends to adsorb to plastic and glass. For these reasons, unused dilute solutions should be discarded and not stored for further use. However, frozen aliquots of a concentrated stock solution (5 µg/µl) in DMSO are expected to be stable for at least a month at -20°C.

4. To the primed RNA from Step 3 in Section 3.1B (6 µl) add the following components and mix by gentle pipetting:

| | |
|--|--------|
| ● (pink) NEBNext First Strand Synthesis Reaction Buffer (5X) | 4 µl |
| ● (pink) Murine RNase Inhibitor | 0.5 µl |
| Actinomycin D (0.1 µg/µl) | 5 µl |
| ● (pink) ProtoScript II Reverse Transcriptase | 1 µl |
| Nuclease free water | 3.5 µl |
| <hr/> | |
| Final volume | 20 µl |

Note: If you are following recommendations in Appendix A, for longer RNA fragments, increase the incubation at 42°C from 15 minutes to 50 minutes in Step 5.

5.  Incubate the sample in a preheated thermal cycler (with the heated lid set at 105°C) as follows:
- 10 minutes at 25°C
 - 15 minutes at 42°C
 - 15 minutes at 70°C
 - Hold at 4°C
6. Proceed directly to Second Strand cDNA Synthesis, Section 3.2.

3.2 Perform Second Strand cDNA Synthesis

- Add the following reagents to the First Strand Synthesis reaction (20 μ l):

| | |
|--|------------|
| Nuclease-free water | 48 μ l |
| ● (orange) Second Strand Synthesis Reaction Buffer (10X) | 8 μ l |
| ● (orange) Second Strand Synthesis Enzyme Mix | 4 μ l |
| Total volume | 80 μ l |
- Mix thoroughly by gentle pipetting.
- Incubate in a thermal cycler for 1 hour at 16°C, with heated lid set at \leq 40°C.

3.3 Purify the Double-stranded cDNA Using 1.8X Agencourt AMPure XP Beads

- Vortex AMPure XP Beads to resuspend.
- Add 144 μ l (1.8X) of resuspended AMPure XP Beads to the second strand synthesis reaction (~80 μ l). Mix well on a vortex mixer or by pipetting up and down at least 10 times.
- Incubate for 5 minutes at room temperature.
- Quickly spin the tube in a microcentrifuge to collect any sample on the sides of the tube. Place the tube on an appropriate magnetic rack to separate beads from supernatant. After the solution is clear (about 5 minutes), carefully remove and discard the supernatant. Be careful not to disturb the beads that contain DNA targets.
- Add 200 μ l of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
- Repeat Step 5 once for a total of 2 washing steps.
- Air dry the beads for 5 minutes while the tube is on the magnetic rack with lid open.

Caution: Do not over dry the beads. This may result in lower recovery of DNA target.

- Remove the tube from the magnet. Elute the DNA target from the beads into 60 μ l 0.1X TE Buffer or 10 mM Tris-HCl. Mix well on a vortex mixer or by pipetting up and down. Quickly spin the tube and incubate for 2 minutes at room temperature. Place the tube in the magnetic rack until the solution is clear.
- Remove 55.5 μ l of the supernatant and transfer to a clean nuclease free PCR tube.



Note: If you need to stop at this point in the protocol samples can be stored at -20°C.

3.4 Perform End Prep of cDNA Library

- Mix the following components in a sterile nuclease free tube:

| | |
|---|--------------|
| Purified double-stranded cDNA (Step 9, Section 2.7) | 55.5 μ l |
| ● (green) NEBNext End Repair Reaction Buffer (10X) | 6.5 μ l |
| ● (green) NEBNext End Prep Enzyme Mix | 3 μ l |
| <hr/> | |
| Total volume | 65 μ l |
- Incubate the sample in a thermal cycler (with the heated lid set at 75°C) as follows:
30 minutes at 20°C
30 minutes at 65°C
Hold at 4°C
- Proceed immediately to Adaptor Ligation.

3.5 Perform Adaptor Ligation

Dilute the ● (red) NEBNext Adaptor* for Illumina (15 μ M) to 1.5 μ M with a 10-fold dilution (1:9) with 10 mM Tris-HCl and 10 mM NaCl for immediate use.

- Add the following components directly to the End Prep Reaction (**Caution: Do not pre-mix the components to prevent adaptor-dimer formation**):

| | |
|------------------------------------|--------------|
| End Prep Reaction | 65 μ l |
| ● (red) Blunt/TA Ligase Master Mix | 15 μ l |
| Diluted NEBNext Adaptor* | 1 μ l |
| Nuclease-free Water | 2.5 μ l |
| <hr/> | |
| Total volume | 83.5 μ l |

*The adaptor is provided in NEBNext Singleplex (NEB #E7350) or NEBNext Multiplex (NEB #E7335, #E7500, #E7710, #E7730, #E6609 or #E7600) Oligos for Illumina.

- Mix by pipetting followed by a quick spin to collect all liquid from the sides of the tube.
- Incubate 15 minutes at 20°C in a thermal cycler.



A precipitate can form upon thawing of the NEBNext Hot Start HiFi PCR Master Mix. To ensure optimal performance, place the master mix at room temperature while purifying the ligation reaction. Once thawed, gently mix by inverting the tube several times.

3.6 Purify the Ligation Reaction Using AMPure XP Beads



Note: If you are selecting for larger size fragments (> 200 nt) follow the size selection recommendations in Appendix A, Chapter 4.

1. Add nuclease-free water to the ligation reaction to bring the reaction volume to 100 μ l. It is important to ensure the final volume is 100 μ l prior to adding AMPure XP Beads.
Note: X refers to the original sample volume of 100 μ l from the above step.
2. Add 100 μ l (1.0X) resuspended AMPure XP Beads and mix well on a vortex mixer or by pipetting up and down at least 10 times.
3. Incubate for 5 minutes at room temperature.
4. Quickly spin the tube in a microcentrifuge and place the tube on an appropriate magnetic rack to separate beads from the supernatant. After the solution is clear (about 5 minutes), discard the supernatant that contain unwanted fragments (**Caution: do not discard the beads**).
5. Add 200 μ l of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
6. Repeat Step 5 once for a total of 2 washing steps.
7. Briefly spin the tube, and put the tube back in the magnetic rack.
8. Completely remove the residual ethanol, and air dry beads for 5 minutes while the tube is on the magnetic rack with the lid open.
Caution: Do not over dry the beads. This may result in lower recovery of DNA target.
9. Remove the tube from the magnet. Elute DNA target from the beads with 52 μ l 0.1X TE or 10 mM Tris-HCl. Mix well on a vortex mixer or by pipetting up and down, incubate for 2 minutes at room temperature. Put the tube in the magnetic rack until the solution is clear.
10. Transfer the 50 μ l supernatant to a clean PCR tube. Discard beads.
11. To the 50 μ l supernatant, add 50 μ l (1.0X) of the resuspended AMPure XP Beads and mix well on a vortex mixer or by pipetting up and down at least 10 times.
12. Incubate for 5 minutes at room temperature.
13. Quickly spin the tube in a microcentrifuge and place the tube on an appropriate magnetic rack to separate beads from the supernatant. After the solution is clear (about 5 minutes), discard the supernatant that contains unwanted fragments (**Caution: do not discard the beads**).

14. Add 200 μ l of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
15. Repeat Step 14 once for a total of 2 washing steps.
16. Briefly spin the tube, and put the tube back in the magnetic rack.
17. Completely remove the residual ethanol, and air dry beads for 5 minutes while the tube is on the magnetic rack with the lid open.

Caution: Do not over dry the beads. This may result in lower recovery of DNA target.

18. Remove the tube from the magnet. Elute DNA target from the beads with 19 μ l 0.1X TE or 10 mM Tris-HCl. Mix well on a vortex mixer or by pipetting up and down, incubate for 2 minutes at room temperature. Put the tube in the magnetic rack until the solution is clear.
19. Without disturbing the bead pellet, transfer 17 μ l of the supernatant to a clean PCR tube and proceed to PCR enrichment.

3.7 PCR Enrichment of Adaptor Ligated DNA



Follow Section 3.7A if you are using the following oligos (10 μ M primer):

NEBNext Singleplex Oligos for Illumina (NEB #E7350)

NEBNext Multiplex Oligos for Illumina (Set 1, NEB #E7335)

NEBNext Multiplex Oligos for Illumina (Set 2, NEB #E7500)

NEBNext Multiplex Oligos for Illumina (Set 3, NEB #E7710)

NEBNext Multiplex Oligos for Illumina (Set 4, NEB #E7730)

NEBNext Multiplex Oligos for Illumina (Dual Index Primers, NEB #E7600)

Follow Section 3.7B if you are using NEBNext Multiplex Oligos for Illumina (96 Index Primers, NEB #E6609).

3.7A PCR Library Enrichment

- To the cDNA (17 μ l) from Step 19 Section 3.6 add the following components and mix by gentle pipetting:

| | |
|---|-------------|
| ● (blue) NEBNext USER Enzyme | 3 μ l |
| ● (blue) NEBNext Q5 Hot Start HiFi PCR Master Mix | 25 μ l |
| ● (blue) Index (X) Primer/i7 Primer*, ** | 2.5 μ l |
| ● (blue) Universal PCR Primer/i5 Primer*, *** | 2.5 μ l |
| <hr/> | |
| Total volume | 50 μ l |

* The primers are provided in NEBNext Singleplex (NEB #E7350) or Multiplex (NEB #E7335, #E7500, #E7710, #E7730, #E7600) Oligos for Illumina. For use with Dual Index Primers (NEB #E7600), look at the NEB #E7600 manual for valid barcode combinations and tips for setting up PCR reactions.

** For use with NEBNext Multiplex Oligos (NEB #E7335, #E7500, #E7710 or #E7730) use only one Index Primer per PCR reaction. For use with Dual Index Primers (NEB #E7600) use only one i7 Primer per reaction.

*** For use with Dual Index Primers (NEB #E7600) use only one i5 Primer per reaction.

- PCR Cycling Conditions

| CYCLE STEP | TEMP | TIME | CYCLES |
|----------------------|------|------------|------------|
| USER Digestion | 37°C | 15 minutes | 1 |
| Initial Denaturation | 98°C | 30 seconds | 1 |
| Denaturation | 98°C | 10 seconds | 12–15*, ** |
| Annealing/Extension | 65°C | 75 seconds | |
| Final Extension | 65°C | 5 minutes | 1 |
| Hold | 4°C | ∞ | |

* The number of PCR cycles should be adjusted based on RNA input. If 10 ng enriched RNA is the starting input, it is recommended to perform 15 cycles of PCR.

** It is important to limit the number of PCR cycles to avoid overamplification.

If overamplification occurs, larger molecular weight products (> 500 bp) will appear on the Bioanalyzer trace.

- Proceed to Section 3.8 (Purify the PCR Reaction using Agencourt AMPure XP Beads).

3.7B PCR Library Enrichment

- To the cDNA (17 μ l) from Step 19 Section 3.6 add the following components and mix by gentle pipetting:

| | |
|---|-----------------------------|
| ● (blue) NEBNext USER Enzyme | 3 μ l |
| ● (blue) NEBNext Q5 Hot Start HiFi PCR Master Mix | 25 μ l |
| ● (blue) Index/ Universal Primer Mix* | 5 μ l |
| Total volume | 50 μl |

- * The primers are provided in NEBNext Multiplex Oligos for Illumina, NEB #E6609. Please refer to the NEB #E6609 manual for valid barcode combinations and tips for setting up PCR reactions.

- PCR Cycling Conditions

| CYCLE STEP | TEMP | TIME | CYCLES |
|----------------------|------|------------|------------|
| USER Digestion | 37°C | 15 minutes | 1 |
| Initial Denaturation | 98°C | 30 seconds | 1 |
| Denaturation | 98°C | 10 seconds | 12–15*, ** |
| Annealing/Extension | 65°C | 75 seconds | |
| Final Extension | 65°C | 5 minutes | 1 |
| Hold | 4°C | ∞ | |

* The number of PCR cycles should be adjusted based on RNA input. If 10 ng enriched RNA is the starting input, it is recommended to perform 15 cycles of PCR.

** It is important to limit the number of PCR cycles to avoid overamplification. If overamplification occurs, larger molecular weight products (> 500 bp) will appear on the Bioanalyzer trace.

- Proceed to Section 3.8 (Purify the PCR Reaction using Agencourt AMPure XP Beads).

3.8 Purify the PCR Reaction using Agencourt AMPure XP Beads

Note: X refers to the original sample volume from the above step.

1. Vortex Agencourt AMPure XP Beads to resuspend.
2. Add 45 μ l (0.9X) of resuspended Agencourt AMPure XP Beads to the PCR reaction (~ 50 μ l). Mix well on a vortex mixer or by pipetting up and down at least 10 times.
3. Incubate for 5 minutes at room temperature.
4. Quickly spin the tube in a microcentrifuge and place the tube on an appropriate magnetic rack to separate beads from the supernatant. After the solution is clear (about 5 minutes), carefully remove and discard the supernatant. Be careful not to disturb the beads that contain DNA targets.
5. Add 200 μ l of freshly prepared 80% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
6. Repeat Step 5 once for a total of 2 washing steps.
7. Air dry the beads for 5 minutes while the tube is on the magnetic rack with the lid open.

Caution: Do not over dry the beads. This may result in lower recovery of DNA target.


8. Remove the tube from the magnetic rack. Elute the DNA target from the beads into 23 μ l 0.1X TE. Mix well on a vortex mixer or by pipetting up and down, quickly spin the tube in a microcentrifuge and incubate for 2 minutes at room temperature. Place the tube in the magnetic rack until the solution is clear.
9. Transfer 20 μ l of the supernatant to a clean PCR tube, and store at -20°C .

3.9 Assess library quality on a Bioanalyzer (Agilent High Sensitivity Chip).

1. Dilute 2–3 μl of the library in 10 mM Tris or 0.1X TE.
2. Run 1 μl in a DNA High Sensitivity chip
3. Check that the electropherogram shows a narrow distribution with a peak size approximately 300 bp.

Note: If a peak at ~ 80 bp (primers) or 128 bp (adaptor-dimer) is shown in the Bioanalyzer traces; Bring up the sample volume (Step 9, Section 3.8) to 50 μl exactly with nuclease-free water and repeat the AMPure XP Bead clean up step (Section 3.8).


Figure 3.1: Example of RNA library size distribution on a Bioanalyzer.



4

Appendix A

NEBNext Ultra Directional
RNA Library Prep Kit for Illumina
Instruction Manual



Note: These recommendations have been optimized using Universal Human Reference Total RNA. Other types of RNA may require different fragmentation times.

Modified fragmentation times for longer RNA inserts.




Figure 4.1: Bioanalyzer traces of RNA as shown in RNA Pico Chip. mRNA isolated from Universal Human Reference RNA (1 μ g) using the NEBNext Poly (A) mRNA Magnetic Isolation Module (NEB #E7490) and Fragmented with First Strand Synthesis Reaction Buffer and Random Primer Mix (2X) at 94°C for 5, 10 or 15 minutes. **For libraries with RNA insert sizes larger than 300 bp, fragment RNA between 5–10 minutes.**

Table 4.1: Recommended size selection conditions for libraries with insert sizes larger than 300 bp.

| LIBRARY PARAMETER | APPROXIMATE INSERT SIZE | 250-400 bp | 300-450 bp | 400-600 bp | 500-700 bp |
|------------------------------------|----------------------------|------------|------------|------------|------------|
| | Approx. Final Library Size | 350-500 bp | 400-550 bp | 500-700 bp | 600-800 bp |
| BEAD VOLUME TO BE ADDED (μ l) | 1st Bead Selection | 45 | 40 | 35 | 30 |
| | 2nd Bead Selection | 20 | 20 | 15 | 15 |

Note: Any differences in insert sizes between the Agilent Bioanalyzer and that obtained from paired end sequencing can be attributed to the higher clustering efficiency of smaller sized fragments.

Size Selection of Adaptor-ligated DNA




For libraries with different size fragment inserts, refer to Table 4.1 for the appropriate volume of beads to be added. The size selection protocol is based on a starting volume of 100 μ l. The protocol below is for libraries with a 300–450 bp insert size.

1. Vortex AMPure XP Beads to resuspend.
2. Adjust the final volume after ligation by adding nuclease free water for a 100 μ l total volume.
3. Add 40 μ l of resuspended AMPure XP Beads to the 100 μ l ligation reaction. Mix well by pipetting up and down at least 10 times.
4. Incubate for 5 minutes at room temperature.
5. Quickly spin the tube and place the tube on an appropriate magnetic stand to separate the beads from the supernatant. After the solution is clear (about 5 minutes), carefully transfer the supernatant containing your DNA to a new tube (**Caution: do not discard the supernatant**). Discard the beads that contain the unwanted large fragments.
6. Add 20 μ l resuspended AMPure XP Beads to the supernatant, mix well and incubate for 5 minutes at room temperature.
7. Quickly spin the tube and place it on an appropriate magnetic stand to separate the beads from the supernatant. After the solution is clear (about 5 minutes), carefully remove and discard the supernatant that contains unwanted DNA. Be careful not to disturb the beads that contain the desired DNA targets (**Caution: do not discard beads**).
8. Add 200 μ l of 80% freshly prepared ethanol to the tube while in the magnetic stand. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
9. Repeat Step 8 once for a total of two washes.
10. Air the dry beads for 5 minutes while the tube is on the magnetic stand with the lid open.

Caution: Do not over dry the beads. This may result in lower recovery of DNA target.

11. Elute the DNA target from the beads into 19 μ l of 10 mM Tris-HCl or 0.1 X TE. Mix well on a vortex mixer or by pipetting up and down. Quickly spin the tube and place it on a magnetic stand. After the solution is clear (about 5 minutes), transfer 17 μ l to a new PCR tube for amplification.

Figure 4.2: Recommended size selection conditions for libraries with insert sizes > 300 bp.



RNA libraries made from Universal Human Reference Total RNA (500 ng) and size selected using different bead/DNA ratios as indicated in Table 4.1. RNA was fragmented at 94°C for 5 minutes.

5

Troubleshooting Guide

NEBNext Ultra Directional
RNA Library Prep Kit for Illumina
Instruction Manual

| OBSERVATIONS | POSSIBLE CAUSES | EFFECT | SUGGESTED SOLUTIONS |
|---|---|---|--|
| Presence of Bioanalyzer peaks <85 bp (Figure 5.1) | <ul style="list-style-type: none"> • Presence of Primers remaining after PCR clean up | Primers cannot cluster or be sequenced, but can bind to flowcell and reduce cluster density | <ul style="list-style-type: none"> • Clean up PCR again with 1.0X AMPure beads (second clean up may result in reduction of library yield) |
| Presence of 127 bp adaptor-dimer Bioanalyzer peak (Figure 5.1) | <ul style="list-style-type: none"> • Addition of non-diluted adaptor • RNA input was too low • RNA was over fragmented or lost during fragmentation • Inefficient Ligation | Adaptor-dimer will cluster and be sequenced. If ratio is low compared to library, may not be a problem but some reads will be dimers. | <ul style="list-style-type: none"> • Dilute adaptor (10 fold dilution) before setting up ligation reaction • Clean up PCR again with 1.0X AMPure beads (second clean up may result in reduction of library yield). |
| Presence of additional Bioanalyzer peak at higher molecular weight than the expected library size (~ 1,000 bp) (Figure 5.2) | <ul style="list-style-type: none"> • PCR artifact (over-amplification). Represents single-stranded library products that have self-annealed. If the PCR cycle number (or PCR input amount) is too high; in the late cycles of PCR the primers become limiting. Therefore, the adaptor sequences on either end of the fragment anneal to each other. This creates molecules with different insert sequences that run slower in the Bioanalyzer. | If ratio is low compared to library, may not be a problem for sequencing | <ul style="list-style-type: none"> • Reduce number of PCR cycles. |
| Broad library size distribution (Figure 5.3) | <ul style="list-style-type: none"> • Under-fragmentation of the RNA | Library size will contain longer insert sizes | <ul style="list-style-type: none"> • Increase RNA fragmentation time |

Figure 5.1:




Figure 5.2:





Figure 5.3:



6

Frequently Asked Questions (FAQs)

NEBNext Ultra Directional
RNA Library Prep Kit for Illumina
Instruction Manual

- Q. What is the difference between the NEBNext Ultra Directional RNA library prep kit for Illumina (E7420) and the NEBNext Ultra RNA library prep kit for Illumina (E7530)?*
- A. The NEBNext Directional RNA library prep workflow preserves information about RNA strand orientation while the NEBNext Ultra RNA library prep does not. The NEBNext Ultra Directional RNA library prep contains dUTP in the second strand synthesis buffer that allows labeling the second strand cDNA and posterior excision with USER enzyme.
- Q. What is the starting material I need to use when preparing libraries using the NEBNext Ultra Directional RNA kit?*
- A. The starting material is Total RNA (100 ng-1 µg); previously isolated mRNA (10-100 ng) or Ribosomal-depleted RNA (10-100 ng).
- Q. Where do I have to start the protocol if I have purified mRNA or Ribosome-depleted RNA?*
- A. If starting material is purified mRNA or ribosome-depleted RNA, proceed to Chapter 3 of the manual.
- Q. Which kit can I use to isolate Poly (A) mRNA from Total RNA?*
- A. To isolate poly (A) mRNA from Total RNA use the NEBNext Poly(A) mRNA Magnetic Isolation Module (NEB #E7490).
- Q. Does the kit provide adaptor and primers?*
- A. No. Adaptors and primers are provided in the NEBNext Singleplex (NEB #E7350) or NEBNext Multiplex (NEB #E7335, #E7500, #E6609, #E7710, #E7730 or #E7600) Oligos for Illumina.

NEBNext First Strand Synthesis Reaction Buffer

#E7421A: 0.192 ml

Concentration: 5X

#E7421AA: 0.768 ml

Store at -20°C

Quality Control Assays

16-Hour Incubation: 50 µl reactions containing 1X First Strand Synthesis Reaction Buffer and 1 µg of HindIII digested Lambda DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 µl reactions containing 1X First Strand Synthesis Reaction Buffer and 1 µg of T3 DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

Endonuclease Activity: Incubation of a 10 µl reaction containing 1X First Strand Synthesis Reaction Buffer with 1 µg of φX174 RF I supercoiled DNA for 4 hours at 37°C results in < 10% conversion to RF II (nicked molecules) as determined by gel electrophoresis.

RNase Activity: Incubation of a 10 µl reaction containing 1X First Strand Synthesis Reaction Buffer with 40 ng of RNA transcript for 16 hours at 37°C resulted in no detectable degradation of RNA as determined by gel electrophoresis.

Phosphatase Activity: Incubation of 1X First Strand Synthesis Reaction Buffer in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl₂) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

Lot Controlled

Random Primers

#E7422A: 0.048 ml

#E7422AA: 0.192 ml

Store at -20°C

Description: This mixture of random hexanucleotides is used to prime DNA synthesis *in vitro* along multiple sites of template RNA.

Sequence: 5' d(N⁶) 3' [N=A,C,G,T]

Phosphorylated: No.

Quality Control Assays

16-Hour Incubation: 50 μ l reactions containing 1 μ l Random Primers and 1 μ g of HindIII digested Lambda DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 μ l reactions containing 1 μ l Random Primers and 1 μ g of T3 DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

Endonuclease Activity: Incubation of a 10 μ l reaction containing 1 μ l Random Primers with 1 μ g of ϕ X174 RF I supercoiled DNA for 4 hours at 37°C results in < 10% conversion to RF II (nicked molecules) as determined by gel electrophoresis.

RNase Activity: Incubation of a 10 μ l reaction containing 1 μ l Random Primers with 40 ng of RNA transcript for 16 hours at 37°C resulted in no detectable RNase activity as determined by gel electrophoresis.

Phosphatase Activity: Incubation of 1 μ l Random Primers in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl₂) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

Lot Controlled

ProtoScript II Reverse Transcriptase

#E7423A: 0.024 ml

Concentration: 200,000 U/ml

#E7423AA: 0.096 ml

Store at -20°C

Description: ProtoScript II Reverse Transcriptase is a recombinant M-MuLV reverse transcriptase with reduced RNase H activity and increased thermostability. It can be used to synthesize first strand cDNA at higher temperatures than the wild type M-MuLV. The enzyme is active up to 50°C, providing higher specificity, higher yield of cDNA and more full-length cDNA product up to 12 kb.

Source: The gene encoding a mutant M-MuLV Reverse Transcriptase (RNase H⁻) is expressed in *E. coli* and purified to near homogeneity.

Supplied in: 20 mM Tris-HCl (pH 7.5), 100 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 0.01% (v/v) IGEPAL® CA-630, 50% (v/v) glycerol

Quality Control Assays

16-Hour Incubation: A 50 µl reaction containing 1 µg of φX174 DNA and 100 units of ProtoScript II Reverse Transcriptase incubated for 16 hours at 37°C resulted in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

Exonuclease Activity: Incubation of a 50 µl reaction containing 100 units of ProtoScript II Reverse Transcriptase with 1 µg of a mixture of single and double-stranded [³H] *E. coli* DNA (10⁵ cpm/µg) for 4 hours at 37°C released < 0.2% of the total radioactivity.

RNase Activity: Incubation of a 10 µl reaction containing 100 units of ProtoScript II Reverse Transcriptase with 40 ng of RNA transcripts for 2 hours at 37°C resulted in no detectable degradation of the RNA as determined by gel electrophoresis.

Phosphatase Activity: Incubation of 100 units of ProtoScript II Reverse Transcriptase in protein phosphatase assay buffer containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

Protein Purity (SDS-PAGE): ProtoScript II Reverse Transcriptase is > 95% pure as determined by SDS PAGE analysis using Coomassie blue detection.

Lot Controlled

Murine RNase Inhibitor

#E7424A: 0.015 ml

#E7424AA: 0.048 ml

Store at -20°C

Description: Murine RNase inhibitor is a 50 kDa recombinant protein of murine origin. The inhibitor specifically inhibits RNases A, B and C. It inhibits RNases by binding noncovalently in a 1:1 ratio with high affinity. It is not effective against RNase 1, RNase T1, S1 Nuclease, RNase H or RNase from *Aspergillus*. In addition, no inhibition of polymerase activity is observed when RNase Inhibitor is used with *Taq* DNA Polymerase, AMV or M-MuLV Reverse Transcriptases, or Phage RNA Polymerases (SP6, T7, or T3).

Recombinant murine RNase inhibitor does not contain the pair of cysteines identified in the human version that is very sensitive to oxidation, which causes inactivation of the inhibitor (1). As a result, murine RNase inhibitor has significantly improved resistance to oxidation compared to the human/porcine RNase inhibitors, even under conditions where the DTT concentration is low. Therefore, it is advantageous to use murine RNase inhibitor in reactions where high concentration DTT is adverse to the reaction (eg. Real-time RT-PCR).

Source: An *E. coli* strain that carries the Ribonuclease Inhibitor gene from mouse.

Supplied in: 20 mM HEPES-KOH, 50 mM KCl, 8 mM DTT and 50% glycerol.

Quality Control Assays

16-Hour Incubation: 50 µl reactions containing a minimum of 40 units of Murine RNase Inhibitor and 1 µg of HindIII digested Lambda DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 µl reactions containing a minimum of 40 units of Murine RNase Inhibitor and 1 µg of T3 DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

Exonuclease Activity: Incubation of a 50 µl reaction containing 200 units of Murine RNase Inhibitor with 1 µg of a mixture of single and double-stranded [³H] *E. coli* DNA (20⁶ cpm/µg) for 4 hours at 37°C released < 0.5% of the total radioactivity.

Latent RNase Assay: Heating the Murine RNase Inhibitor for 20 minutes at 65°C, followed by incubation of a 10 µl reaction containing 40 units of RNase Inhibitor with 40 ng of RNA transcript for 4 hours at 37°C resulted in no detectable degradation of RNA as determined by gel electrophoresis.

RNase Activity: Incubation of a 10 µl reaction containing 40 units of Murine RNase Inhibitor with 40 ng of RNA transcript for 16 hours at 37°C resulted in no detectable degradation of RNA as determined by gel electrophoresis.

Endonuclease Activity: Incubation of a 10 µl reaction containing 40 units of Murine RNase Inhibitor with 1 µg of φX174 RF I supercoiled DNA for 4 hours at 37°C results in < 10% conversion to RF II (nicked molecules) as determined by gel electrophoresis.

Phosphatase Activity: Incubation of a minimum of 40 units of Murine RNase Inhibitor in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl₂) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

Lot Controlled

References:

1. Kim, B.M. et al. (1999). *Protein Science*, 8, 430-434.

NEBNext Second Strand Synthesis Enzyme Mix

#E7425A: 0.096 ml

#E7425AA: 0.384 ml

Store at -20°C

Description: NEBNext Second Strand Synthesis Enzyme Mix is optimized to convert short single-stranded cDNAs to double-stranded cDNAs.

Quality Control Assays

SDS-PAGE Purity: SDS-PAGE analysis of this enzyme indicates > 95% enzyme purity.

Endonuclease Activity: Incubation of a 10 μl reaction containing 1 μl Second Strand Synthesis Enzyme Mix with 1 μg of ϕX174 RF I supercoiled DNA for 4 hours at 37°C results in < 10% conversion to RF II (nicked molecules) as determined by gel electrophoresis.

Phosphatase Activity: Incubation of a minimum of 1 μl Second Strand Synthesis Enzyme Mix in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl_2) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

Functional Activity: One unit of the *E. coli* DNA Ligase ligated 50% of HindIII fragments of λ DNA (5' DNA termini concentration of 0.12 μM , 300 $\mu\text{g}/\text{ml}$) in a total reaction volume of 20 μl in 30 minutes at 16°C in 1X *E. coli* DNA Ligase Reaction Buffer. One unit of *E. coli* DNA Polymerase I incorporated 10 nmol of dNTP into acid-insoluble material in a total reaction volume of 50 μl in 30 minutes at 37°C in 1X EcoPol Reaction Buffer with 33 μM dNTPs including [^3H]-dTTP and 70 $\mu\text{g}/\text{ml}$ denatured herring sperm DNA. Incubation of 50 units of RNase H with 1 μg sonicated and denatured [^3H]-DNA (10^5 cpm/ μg) for 30 minutes at 37°C in 50 μl reaction buffer released < 0.1% radioactivity.

Lot Controlled

Reference:

1. Gubler et al. (1983). *Gene* 25, 263–269.

NEBNext Second Strand Synthesis Reaction Buffer with dUTP Mix

#E7426A: 0.192 ml

Concentration: 10X

#E7426AA: 0.768 ml

Store at -20°C

Quality Control Assays

16-Hour Incubation: 50 μl reactions containing 1X NEBNext Second Strand Synthesis Reaction Buffer with dUTP mix and 1 μg of HindIII digested Lambda DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 μl reactions containing 1X Second Strand Synthesis Reaction Buffer and 1 μg of T3 DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

Endonuclease Activity: Incubation of a 10 μl reaction containing 1X NEBNext Second Strand Synthesis Reaction Buffer with dUTP mix with 1 μg of ϕX174 RF I supercoiled DNA for 4 hours at 37°C results in $< 10\%$ conversion to RF II (nicked molecules) as determined by gel electrophoresis.

RNase Activity: Incubation of a 10 μl reaction containing 1X NEBNext Second Strand Synthesis Reaction Buffer with dUTP mix with 40 ng of RNA transcript for 16 hours at 37°C resulted in no detectable degradation of RNA as determined by gel electrophoresis.

Phosphatase Activity: Incubation of 1X NEBNext Second Strand Synthesis Reaction Buffer with dUTP mix in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl_2) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

Lot Controlled

Limited Use Label License: This product and its use are subject to one or more issued and/or pending U.S. and foreign patent applications owned by Max Planck Gesellschaft. The purchase of this product from New England Biolabs, Inc., its affiliates, or its authorized resellers and distributors conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for profit entity) or for limited commercial purposes that are pre-approved by New England Biolabs, Inc. The purpose of this product does not convey a license under any claims in the foregoing patents or patent applications directed to producing the product. The buyer cannot sell or otherwise transfer this product or its components to a third party or otherwise use this product for the following COMMERCIAL PURPOSES: (1) use of the product or its components in manufacturing; or (2) use of the product or its components for therapeutic or prophylactic purposes in humans or animals.

NEBNext End Prep Enzyme Mix

#E7371A: 0.072 ml

#E7371AA: 0.288 ml



Store at -20°C

Description: NEBNext End Prep Enzyme Mix is optimized to convert 5 ng–1 µg of fragmented DNA to repaired DNA having 5´-phosphorylated dA-tailed ends.

Quality Control Assays

SDS-PAGE Purity: SDS-PAGE analyses of each individual enzyme indicates > 95% enzyme purity.

Endonuclease Activity: Incubation of a minimum of 10 µl of this enzyme mix with 1 µg of φX174 RF I DNA in assay buffer for 4 hours at 37°C in 50 µl reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

Phosphatase Activity: Incubation of a minimum of 10 µl of this enzyme mix in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl₂) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

Functional Activity (Nucleotide Incorporation, Phosphorylation and dA-Tailing): 1 µl of this enzyme mix repairs and phosphorylates the ends of > 95% of 0.5µg of DNA fragments containing both 3´ and 5´ overhangs with 20 minutes at 25°C, in 1X End Repair Reaction buffer, as determined by capillary electrophoresis.

Lot Controlled

NEBNext End Repair Reaction Buffer

#E7372A: 0.156 ml

Concentration: 10X

#E7372AA: 0.624 ml

Store at -20°C

Quality Control Assays

16-Hour Incubation: 50 µl reactions containing this reaction buffer at 1X concentration and 1 µg of HindIII digested Lambda DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 µl reactions containing this reaction buffer at 1X concentration and 1 µg T3 DNA incubated for 16 hours at 37°C also results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

Endonuclease Activity: Incubation of this reaction buffer at a 1X concentration with 1 µg of φX174 RF I DNA for 4 hours at 37°C in 50 µl reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

RNase Activity: Incubation of this reaction buffer at 1X concentration with 40 ng of a FAM-labeled RNA transcript for 16 hours at 37°C results in no detectable RNase activity as determined by polyacrylamide gel electrophoresis.

Phosphatase Activity: Incubation of this reaction buffer at a 1X concentration in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl₂) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

Lot Controlled

Blunt/TA Ligase Master Mix

#E7373A: 0.360 ml

#E7373AA: 0.720 ml (2 vials provided)

Store at -20°C

Description: Blunt/TA Ligase Master Mix is a ready-to-use solution of T4 DNA Ligase, proprietary ligation enhancer, and optimized reaction buffer.

Quality Control Assays

16-Hour Incubation: 50 µl reactions containing this reaction buffer at 1X concentration and 1 µg of HindIII digested Lambda DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 µl reactions containing this reaction buffer at 1X concentration and 1 µg T3 DNA incubated for 16 hours at 37°C also results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

Endonuclease Activity: Incubation of this reaction buffer at a 1X concentration with 1 µg of φX174 RF I DNA for 4 hours at 37°C in 50 µl reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

RNase Activity: Incubation of this reaction buffer at 1X concentration with 40 ng of a FAM-labeled RNA transcript for 16 hours at 37°C results in no detectable RNase activity as determined by polyacrylamide gel electrophoresis.

Transformation Assay: LITMUS™ 28 vector is cut with EcoRV (blunt), treated with calf intestinal phosphatase and gel purified. Blunt inserts from a HaeIII digest of φX174 DNA are ligated into the vector at a 3:1 insert:vector ratio using the Blunt/TA Ligase Master Mix Protocol. Ligation products are transformed as described.

Each lot exceeds the following standards:

| | Efficiency (transformants/µg) | |
|--------------|--------------------------------------|-------------------------|
| | Recircularization | Insertion |
| Blunt ends | > 1 x 10 ⁷ | > 2.5 x 10 ⁶ |
| Uncut vector | > 1 x 10 ⁸ | |

Lot Controlled

Nuclease-free Water

#E7431A: 8 ml

#E7431AA: 30 ml

Store at -20°C or 4°C

Description: Nuclease-free Water is free of detectable DNA and RNA nucleases and phosphatases and suitable for use in DNA and RNA applications.

Quality Control Assays

16-Hour Incubation: 50 μl reactions containing Nuclease-free Water and 1 μg of HindIII digested Lambda DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 μl reactions containing Nuclease-free Water and 1 μg of T3 DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

Endonuclease Activity: Incubation of a 10 μl reaction containing Nuclease-free Water with 1 μg of ϕX174 RF I supercoiled DNA for 4 hours at 37°C results in $< 10\%$ conversion to RF II (nicked molecules) as determined by gel electrophoresis.

RNase Activity: Incubation of a 10 μl reaction containing Nuclease-free Water with 40 ng of RNA transcript for 16 hours at 37°C resulted in no detectable degradation of RNA as determined by gel electrophoresis.

Phosphatase Activity: Incubation of 1X Second Strand Synthesis Reaction Buffer in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl_2) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

Lot Controlled

NEBNext Q5 Hot Start HiFi PCR Master Mix

E6625A: 0.6 ml

Concentration: 2X

E6625AA: 1.2 ml (2 vials provided)

Store at -20°C

Description: The NEBNext Q5 Hot Start HiFi PCR Master Mix is specifically optimized for robust, high-fidelity amplification of next-generation sequencing (NGS) libraries, regardless of GC content. The polymerase component of the master mix, Q5 High-Fidelity DNA Polymerase, is a novel thermostable DNA polymerase that possesses 3'→5' exonuclease activity, and is fused to a processivity-enhancing Sso7d domain. Q5 also has an ultra-low error rate (> 100-fold lower than that of *Taq* DNA Polymerase and ~12-fold lower than that of *Pyrococcus furiosus* (Pfu) DNA Polymerase). The buffer component of the master mix has been optimized for robust amplification, even with GC-rich amplicons and offers enhanced compatibility with a variety of beads used in typical NGS workflows. These features make the NEBNext Q5 Hot Start HiFi PCR Master Mix ideal for NGS library construction. This convenient 2X master mix contains dNTPs, Mg⁺⁺ and a proprietary buffer, and requires only the addition of primers and DNA template for robust amplification. The inclusion of the hot start aptamer allows convenient room temperature reaction set up.

Quality Control Assays

16-Hour Incubation: A 50 µl reaction containing NEBNext Q5 Hot Start HiFi PCR Master Mix and 1 µg of HindIII digested λ DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 µl reactions containing NEBNext Q5 Hot Start HiFi PCR Master Mix and 1 µg of T3 DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

Phosphatase Activity: Incubation of NEBNext Q5 Hot Start HiFi PCR Master Mix in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl₂) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

Functional Activity (Multiplex PCR, Bead Inhibition): 30 cycles of PCR amplification of 20 ng genomic DNA with and without carboxylated magnetic beads in a 50 µl reaction containing 0.5 µM 4-plex primer mix and 1X NEBNext Q5 Hot Start HiFi PCR Master Mix result in the four expected amplicons and no inhibition of amplification in the presence of the beads.

Lot Controlled

This product is covered by one or more Patents.

This product is licensed from Bio-Rad Laboratories, Inc. under U.S. Pat. Nos. 6,627,424; 7,541,170; 7,670,808; 7,666,645 and corresponding patents in other countries for use only in: (a) standard (non-real time) PCR in the research field only, but not real-time PCR or digital PCR; (b) any *in vitro* diagnostics application, except for applications using real-time or digital PCR; and (c) any non-PCR applications in DNA sequencing, isothermal amplification and the production of synthetic DNA.

NEBNext USER Enzyme

#E7428A: 0.072 ml

#E7428AA: 0.288 ml

Store at -20°C

Supplied in: 50 mM KCl, 5 mM NaCl, 10 mM Tris-HCl (pH 7.4 @ 25°C), 0.1 mM EDTA, 1 mM DTT, 175 µg/ml BSA and 50% Glycerol.

Quality Control Assays

Non-Specific DNase Activity (16 Hour): A 50 µl reaction in NEBuffer 1 containing 1 µg of Lambda DNA and a minimum of 50 units of Uracil DNA Glycosylase incubated for 16 hours at 37°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis. A 50 µl reaction in Endonuclease VIII Reaction Buffer containing 1 µg of Lambda-HindIII DNA and a minimum of 25 units of Endonuclease VIII incubated for 16 hours at 37°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

Exonuclease Activity (Radioactivity Release): A 50 µl reaction in NEBuffer 1 containing 1 µg of a mixture of single and double-stranded [³H] *E. coli* DNA and a minimum of 50 units of Uracil DNA Glycosylase incubated for 4 hours at 37°C releases < 0.1% of the total radioactivity. A 50 µl reaction in Endonuclease VIII Reaction Buffer containing 1 µg of a mixture of single and double-stranded [³H] *E. coli* DNA and a minimum of 10 units of Endonuclease VIII incubated for 4 hours at 37°C releases < 0.5% of the total radioactivity.

Endonuclease Activity (Nicking): A 50 µl reaction in UDG Reaction Buffer containing 1 µg of supercoiled φX174 DNA and a minimum of 50 units of Uracil DNA Glycosylase incubated for 4 hours at 37°C results in < 10% conversion to the nicked form as determined by agarose gel electrophoresis.

Phosphatase Activity: Incubation of a minimum of 10 µl of USER at a 1X concentration in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl₂) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

Lot Controlled

Revision History:

| Revision # | Description |
|------------|--|
| 2.0 | Added RNA input recommendations, removed the size selection for 200 bp fragments - replaced with clean up step. Added additional recommendation for larger insert sizes (Appendix A), Troubleshooting Guide, FAQs. Removed additional washing step in PolyA Isolation Protocol. Moved stopping point from after Second Strand cDNA Synthesis to follow the clean up step.Changed First Strand cDNA Synthesis conditions from 50 minutes at 42°C to 15 minutes at 42°C. Added recommendation to dilute the NEBNext adaptor. |
| 2.1 | Renamed "elution buffer" in text to "Tris Buffer". |
| 3.0 | Added protocol for use with NEB #E6310. Changed NEBNext adaptor dilution recommendations. Changed AMPure Bead drying time to 5 minutes. Updated final library elution and dilution for Bioanalyzer to 10 mM Tris or 0.1X TE. Changed ratio of AMPure Beads to 0.9X in final clean up after PCR reaction. |
| 4.0 | Included protocol for use with NEBNext Q5 Hot Start HiFi PCR Master Mix. Include protocol for changes in concentration of NEBNext Singleplex and Multiplex Oligos for Illumina. Changed all AMPure Bead elutions to 0.1X TE or 10 mM Tris-HCl. Added 2 minute incubation after eluting DNA from AMPure Beads. |
| 5.0 | Removed protocol for use with NEBNext High-Fidelity 2X PCR Master Mix. Include protocol for use with NEBNext Multiplex Oligos (96 Index Primers, NEB #E6609). |
| 6.0 | Volume of beads increased from 15 µl to 20 µl. Additional mixing and incubation steps were added after each thermocycler incubation. Elution volume changed before fragmentation. Updated table on page 40. |
| 7.0 | Component change: The name, part # and formulation of RNase H has changed. |
| 8.0 | Protocol updated to include NEB #E7710 and NEB #E7730. Correction of typos and clarifications in several places. Section C in the PCR setup step was removed because all of the 25 µM primers are now expired. |



USA

New England Biolabs, Inc.
240 County Road
Ipswich, MA 01938-2723
Telephone: (978) 927-5054
Toll Free: (USA Orders) 1-800-632-5227
Toll Free: (USA Tech) 1-800-632-7799
Fax: (978) 921-1350
e-mail: info@neb.com
www.neb.com

CANADA

New England Biolabs, Ltd.
Telephone: (905) 665-4632
Toll Free: 1-800-387-1095
Fax: (905) 665-4635
Fax Toll Free: 1-800-563-3789
e-mail: info.ca@neb.com
www.neb.ca

CHINA, PEOPLE'S REPUBLIC

New England Biolabs (Beijing), Ltd.
Telephone: 010-82378265/82378266
Fax: 010-82378262
e-mail: info@neb-china.com
www.neb-china.com

FRANCE

New England Biolabs France
Free Call: 0800-100-632
Free Fax: 0800-100-610
e-mail: info.fr@neb.com
www.neb-online.fr

GERMANY & AUSTRIA

New England Biolabs GmbH
Telephone: +49/(0)69/305 23140
Free Call: 0800/246 5227 (Germany)
Free Call: 00800/246 52277 (Austria)
Fax: +49/(0)69/305 23149
Free Fax: 0800/246 5229 (Germany)
e-mail: info.de@neb.com
www.neb-online.de

JAPAN

New England Biolabs Japan, Inc.
Telephone: +81 (0)3 5669 6191
Fax: +81 (0)3 5669 6192
e-mail: info.jp@neb.com
www.nebj.jp

SINGAPORE

New England Biolabs Pte. Ltd.
Telephone: +65 6776 0903
Fax: +65 6778 9228
e-mail: sales.sg@neb.com
www.neb.sg

UNITED KINGDOM

New England Biolabs (UK) Ltd.
Telephone: (01462) 420616
Call Free: 0800 318486
Fax: (01462) 421057
Fax Free: 0800 435682
e-mail: info.uk@neb.com
www.neb.uk.com